CATALOGUE Ver. 13

DDS

Drug Delivery Systems

Table of Contents

| | |
|---------------------------------------|---|
| NOF CORPORATION OVERVIEW | 3 |
| DDS PRODUCT OVERVIEW | 4 |

I . Activated PEG for PEGylation 6

NOF's Product Advantages

- High-purity mPEG-OH, Starting Material for Activated PEGs
- About PEGylation
- Capabilities

High-purity Activated PEGs

- Application Chart for PEG Derivatives
- Custom Manufacturing
- 1. Monofunctional Linear PEGs
- 2. Bi-Functional PEGs
- 3. Multi-Arm PEGs
 - 3-1. 4-arm Functional PEG Series
 - 3-2. 8-arm Functional PEG Series
- 4. Branched PEGs
 - 4-1. 2-arm Branched PEG
 - 4-2. 3-arm Branched PEG
 - 4-3. 4-arm Branched PEG
 - 4-4. Lysine Branched PEG
- 5. Heterofunctional PEGs
- 6. Forked PEGs
- 7. Releasable PEGs
- 8. Comb-shaped co-polymers

II . Phospholipids and Lipids for Liposomal Formulations 28

1. Purified Phospholipids from natural sources
2. Highly-purified Synthetic Phospholipids
3. PEGylated Lipids
4. Functionalized Phospholipids
5. Novel Lipids and Cationic Lipids
6. Empty Liposomes <COATSOME® EL Series>
7. Transferrin Pre-liposome
8. Liposomal Formulations
9. LIPONIZER™ for Liposome Production
10. NOF's Capabilities and Customer Advantages

III . Ultra-purity Polysorbate 80 51

IV . Novel Excipients for Pharmaceutical Development

| | |
|--|----|
| PUREBRIGHT® MB Series | 53 |
| Self-emulsifying Solubilizer SL-11 | 55 |
| PUREBRIGHT® SL Series | 57 |
| NOFABLE™ Series | 59 |
| SUNBRIGHT® DKH-02HB, DKH-03HB and DKH-04HB | 60 |
| SUNBRIGHT® OE Series | 61 |

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The NOF Group, pursues multi-faceted business development in nine divisions of activities based on its own technologies, NOF endeavors to realize its management philosophy of “creating new values in a broad spectrum from biosphere to outer space” by focusing particularly on the fields of “life science”, “electronics and information” and “environment and energy”

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Sales : 1.8 billion US\$ (Fiscal year 2010)
US\$=¥86.81

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

● Corporate Organization

DDS Development Division is organized within NOF's family of complementary business units, offering innovation and cGMP capabilities to pharmaceutical and biomedical products.

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● Overseas Subsidiaries and Joint Ventures

DDS PRODUCT OVERVIEW

I . Activated PEG for PEGylation

We manufacture high-purity methoxy-polyethylene glycols (mPEGs) with molecular weights from 2kDa to 80kDa, producing high-purity activated PEG derivatives of methoxy PEG amines, maleimides and carboxylic acids. Various alterations of the terminal PEG chemistry permit drug modification as per customers' expectations and needs. Furthermore, use of PEGylated phospholipids composed of our high-purity phospholipids and PEGs prolong the plasma circulating half-life of liposome drugs.

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II . Phospholipids and Lipids for Liposomal Formulations

In early 1990's, we commenced supply of GMP-grade phospholipids to every corner of the world. Our products have been used in commercialized liposomal drugs for anti-fungal and oncology use, which led to our outstanding reputation in this field. For the benefit of our customers, we have also developed a ready-to-use liposomal formulation, called EMPTY LIPOSOMES, into which customers can encapsulate their own drugs.

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III . Ultra-purity Polysorbate 80

Our Polysorbate 80(HX2) is the ultra-purity polysorbate 80 surfactant, and we are proud of its rating of the highest quality around the world, with less impurities compared to conventional grade. Favorable features, such as less impurities and lower susceptibility to oxidation, are responsible for negligible production of peroxides and aldehydes, which are well known to induce degradation of drug formulations. In addition, owing to its high quality, Polysorbate 80(HX2) shows superior safety features, including a lower incidence of allergic, hemolytic, cytotoxic and acute toxicity reactions compared to conventional grade.

Polysorbate 80(HX2) can be applied to various drugs as a multi-purpose surfactant, e.g., as a solubilizing agent for poorly hydrophilic drugs, as a stabilizing agent for injectable preparations containing proteins, as a stabilizer for vaccines.

Polysorbate 80(HX2) is vegetable-derived material and meets the ICH Harmonized Tripartite Guideline on New Drug Substances and Products (published, April 2006). This has validated its extensive use around the world as a standard product in the pharmaceutical field.

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IV . Novel Excipients for Pharmaceutical Development

PUREBRIGHT[®] MB series

MPC[™] Polymers for Hydrophobic Drug Formulations

Through introduction of a hydrophilic phosphorylcholine moiety into a hydrophobic methacrylic acid polymer, we offer innovative formulation excipients for any hydrophobic drugs. Using our PUREBRIGHT MB series, any drug can be dissolved in aqueous solutions.

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Self-emulsifying Solubilizer SL-11 Novel Nano Formulation for Hydrophobic Drugs

Solubilizer SL-11 is clear liquid with superior solvent properties for many kinds of hydrophobic drugs to prepare nano emulsions with particle sizes below 50 nm. Solubilizer SL-11 is highly suitable as a SEDDS (Self Emulsifying Drug Delivery System) agent. Solubilizer SL-11 can also help increase intestinal absorption. Nano emulsions and SEDDS agents containing hydrophobic drugs can be easily prepared following simple protocols.

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PUREBRIGHT® SL Series Hydrophobic Drug Solubilization Kit

We supply new solubilization kits for hydrophobic drugs by using our proprietary DDS technologies. Because of the simple procedures for kit use, screening of new chemical entities and safety studies can be easily performed.

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NOFABLE™ Series Ultra-Purity Oleic Acid and Derivatives

Using high-purity oleic acid, we produce the non-ionic surfactant (NOFABLE™ series) and its derivatives. Since impurities derived from oxidation are eliminated in our process for the product, any degradation and degeneration of drug formulations using it can be prevented. In addition, oil adjuvant formulation for vaccines also exploits the NOFABLE™ series.

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SUNBRIGHT® DKH-02HB, DKH-03HB and DKH-04HB MACROGOL (PEG-200, 300 and 400)

We have been providing ultra-pure polyethylene glycol (PEG 200, 300 and 400) as an excipient specifically useful for pharmaceutical formulations without any impurities such as EG (ethylene glycol) and DEG (diethylene glycol). We supply PEG200, 300 and 400 that meet the requirements stipulated in the USP, EP and JP /JPE.

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SUNBRIGHT® OE Series Biocompatible PEG Anchors for Cell Membrane Insertion [BAM]


We provide a cell-modifying agent that acts without exerting any injurious effects on the cells. Employment of the BAM concept contributes to both modification of the cells or tissue surfaces with physiologically active substances, such as proteins or drugs, and live cell immobilization of cells or tissues on the surfaces of various kinds of materials without any damage to the target cells and tissues.

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Using these high-quality products, NOF CORPORATION provides custom formulation expertise and assistance on request.

I

Activated PEG for PEGylation <SUNBRIGHT® Series>

NOF has a long reputation for producing a high-quality Activated PEGs -the SUNBRIGHT Series- that possess the most suitable quality for preparation of physiologically active drugs and biologics. Our SUNBRIGHT series is characterized as extremely high-purity PEGs with a variety of functional groups. The products of this series are internationally recognized as optimal PEGylation drugs.

As illustrated in this figure, application of Activated PEGs with various functional groups enables introduction of PEG chains into drugs, enzymes, phospholipids and other biologics. Covalent conjugation of hydrophobic macromolecules with Activated PEGs leads to the formation of macromolecular micelles (Polymer Micelles), which allow homogeneous dispersion of hydrophobic drugs in aqueous media. Furthermore, when PEGylated phospholipids are used as liposomes, the aqueous corona added stabilizes homogeneous dispersal of the liposome preparations encapsulating drugs within them.

PEGylation plays an important role in the stabilization of drugs increased circulation time, reduction of their antigenicity and decrease in drug dosing, besides augmenting the targeting ability via binding biologics onto their surfaces.

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
NOF's Product Advantages

High-purity mPEG-OH, Starting Material for Activated PEGs

NOF manufactures high-purity methoxy-PEG-OH characterized by absence of contaminating by impurities and narrow molecular weights distribution. These highly pure methoxy-PEG-OH have so far been used for various commercial PEGylation drugs, such as PEG-Interferon and PEG-GCSF.

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

NOF CORPORATION is the sole worldwide patent holder of these high-purity mPEG-OH and all their derivatives.
※ Patent No. : JP3050228, U.S. Patent 6455639, KR 0358276

Our methoxy-PEG-OH have the advantage of extremely low diol contents relative to those of competitors' products. As shown in the following graphs, our methoxy-PEG-5000 contains remarkably low levels of diol as impurities, besides providing a narrower distribution of molecular weights in comparison with our competitors' products. Employment of our methoxy-PEG-OH as the starting materials yields higher purity and higher activity of Activated PEGs. We manufacture highly pure methoxy-PEGs with molecular weights from 2kDa to 80kDa.

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Comparisons of GPC analysis of NOF mPEG-OH (MW5,000) (left) and conventional product (right).

About PEGylation

PEGylation technology is applied for biologics modifications including cytokines, therapeutic antibody fragments, interleukins, hormones, oligonucleotides and other proteins and peptides.

PEGylation can be effective for :

- Improving Bioavailability
- Prolonging Blood Circulation
- Maximizing Pharmacokinetics
- Low-profile Immunogenicity

NOF provides high-performance activated PEGs for PEGylated drugs from early development stage to commercial use with our 20 years of experience of manufacturing high-quality methoxy polyethylene glycol. Our starting material for activated PEGs for pharmaceutical area. Our activated PEGs have a narrow polydispersity and low diol content with wide range of molecular weights, from 2kDa to 80kDa.

PEGylation Service

Please feel free to contact us if you need PEGylated drugs using NOF Activated PEGs, SUNBRIGHT® series. Contact : ddsinfo@nof.co.jp

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
Reference

1) PEG-conjugation

Following figure shows several examples of NOF PEG-aldehyde conjugation with human insulin for various PEG molecular weights.

Analysis of PEG-ALD-Insulin by SDS-PAGE with Coomassie brilliant blue. Insulin was reacted with NOF PEG-ALD at 4°C. The PEGylated insulin was purified on Q Sepharose Fast Flowcolumn.

lane M molecular weight markers
lane 1 insulin
lane 2 PEG(ME-050AL)-insulin eluate
lane 3 PEG(ME-100AL)-insulin eluate
lane 4 PEG(ME-200AL)-insulin eluate
lane 5 PEG(ME-300AL)-insulin eluate

2) Coupling of PEG derivatives to Proteins: examples of simple conjugation conditions

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1. Coupling of PEG-NHS derivatives to protein amines (PEG-NHS + Protein-NH₂)
Conditions 1: 50 mM phosphate buffer (pH 7.2), 4°C, 6 hrs^{*1)}
Conditions 2: Borate-phosphate buffer (pH 8.0), 25°C, 2 hrs^{*2)}
2. Coupling of PEG-Aldehyde derivatives to protein amines (PEG-Ald + Protein-NH₂)
Conditions : sodium cyanoborohydride (10 equiv.), 4°C, 20 hrs^{*3)}
3. Coupling of PEG-Maleimide derivatives to protein sulfhydryls (PEG-Maleimide + Protein-SH)
Conditions: 100 mM phosphate buffer (pH 6.5), 4°C, 4 hrs
4. Coupling of PEG-NH₂ derivatives to protein carboxylates (PEG-NH₂ + Protein-COOH)
Conditions: 50mM phosphate buffer (pH 7.2), WSC (2 equiv.), 4°C, 10 hrs
5. Coupling of PEG-p-Nitrophenyloxycarbonyl derivatives to protein amines (PEG-NP + Protein-NH₂)
Conditions: borate-phosphate buffer (pH 8.0-8.3), r.t, overnight^{*2)}

References: *1) Abuchowski A. et al., *Cancer Biochem. Biophys.* **7**,175(1984)
*2) Sartore L. et al., *Appl.Biochem. Biotechnol.* **27**,45(1991)
*3) U.S.Patent 5,824,784

Capabilities

cGMP Manufacturing Facilities

NOF Activated PEGs, SUNBRIGHT® series are produced in facilities, using state-of-the-art technology, operated under cGMP. The cGMP facilities are capable of producing small to large batch sizes, using proprietary know-how with scale-up procedures, depending on customers' needs.

These cGMP facilities have been audited by pharmaceutical companies throughout the world. NOF receives a consistent high reputation from our customers.

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Research & Development

Our new R&D facility was opened in December 2005. This allows us to continue our development of novel activated PEGs and new technologies for PEGylation.

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Analytical Technologies

NOF has more than 20 years history for manufacturing of high-quality mPEGs and Activated PEGs. These experiences endow NOF with a thorough knowledge of PEGs and considerable achievement in analysis of Activated PEGs. We also invented innovative analytical methods for assaying purity of activated PEGs and for determining impurity levels.

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High-purity Activated PEG

Application Chart for PEG Derivatives

We supply different series of PEG derivatives with various versatile functional groups.

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* PEG Reactive group (X)
-O-(CH₂CH₂O)_n-X or -O-X

Custom Manufacturing

NOF has over 20 years of experience in supplying high-quality methoxy polyethylene glycol (mPEG) to pharmaceutical markets. The highest quality mPEG is our starting material of PEG derivatives, patented through our proprietary technology.

Based on our extensive range of technological know-how and capabilities, NOF offers custom manufacturing from early stage development to commercial scale.

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1. Monofunctional Linear PEGs

The mono-functional SUNBRIGHT® series comprising highly purified methoxy PEG as the starting material contains negligible bifunctional PEG derived from diol impurities. Accordingly, bridging reactions (sidereaction), an obstacle to PEGylation can be completely prevented.

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
•SUNBRIGHT® CS, GS, AS, HS and TS Series (NHS active esters/carbonate)

We supply five types of monomethoxy-NHS activated ester/carbonate-PEGs that differ in terms of the NHS carboxylate attachment chemistry. The reactivity order of these linkers in aqueous solution from lowest to highest is HS<TS<GS<CS<AS. After PEGylation, CS and GS types possessing ester groups are relatively susceptible to hydrolysis, whereas AS, TS and HS types are resistant to hydrolysis because of their different chemical structures. Therefore, our customers can select the appropriate type according to the purposes among the five types available.

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
* : make-to-order


GS Type

| Product name | MW |
|--------------------|--------|
| SUNBRIGHT ME-050GS | 5,000 |
| SUNBRIGHT ME-100GS | 10,000 |
| SUNBRIGHT ME-200GS | 20,000 |
| SUNBRIGHT ME-300GS | 30,000 |
| SUNBRIGHT ME-400GS | 40,000 |


New

AS Type


| Product name | MW |
|--------------------|-------|
| SUNBRIGHT ME-020AS | 2,000 |
| SUNBRIGHT ME-050AS | 5,000 |

HS Type


| Product name | MW | |
|--------------------|--------|-----|
| SUNBRIGHT ME-050HS | 5,000 | |
| SUNBRIGHT ME-100HS | 10,000 | New |
| SUNBRIGHT ME-200HS | 20,000 | |
| SUNBRIGHT ME-300HS | 30,000 | |
| SUNBRIGHT ME-400HS | 40,000 | |

TS Type


| Product name | MW | |
|--------------------|--------|--|
| SUNBRIGHT ME-050TS | 5,000 | |
| SUNBRIGHT ME-100TS | 10,000 | |
| SUNBRIGHT ME-120TS | 12,000 | |
| SUNBRIGHT ME-200TS | 20,000 | |
| SUNBRIGHT ME-300TS | 30,000 | |
| SUNBRIGHT ME-400TS | 40,000 | |

•SUNBRIGHT® NP Series (p-Nitrophenyl Carbonate PEG)

| Product name | MW | |
|--------------------|--------|--|
| SUNBRIGHT MENP-20H | 2,000 | |
| SUNBRIGHT MENP-50H | 5,000 | |
| SUNBRIGHT MENP-10T | 10,000 | |
| SUNBRIGHT MENP-20T | 20,000 | |
| SUNBRIGHT MENP-30T | 30,000 | |
| SUNBRIGHT MENP-40T | 40,000 | |

•SUNBRIGHT® AL Series (Aldehyde PEG)


Exploiting our latest technologies, we supply mPEG-aldehydes without any reactive problematic impurities (e.g., bifunctional aldehydes).

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| Product name | MW | |
|---------------------|--------|--|
| SUNBRIGHT ME-050AL | 5,000 | |
| SUNBRIGHT ME-100AL | 10,000 | |
| SUNBRIGHT ME-200AL | 20,000 | |
| SUNBRIGHT ME-300AL | 30,000 | |
| SUNBRIGHT ME-400AL2 | 40,000 | |


•SUNBRIGHT® PA Series (Aminopropyl PEG)

| Product name | MW |
|--------------------|--------|
| SUNBRIGHT MEPA-20H | 2,000 |
| SUNBRIGHT MEPA-50H | 5,000 |
| SUNBRIGHT MEPA-10T | 10,000 |
| SUNBRIGHT MEPA-12T | 12,000 |
| SUNBRIGHT MEPA-20T | 20,000 |
| SUNBRIGHT MEPA-30T | 30,000 |
| SUNBRIGHT MEPA-40T | 40,000 |

New


•SUNBRIGHT® EA Series (Aminoethyl PEG)

| Product name | MW |
|--------------------|--------|
| SUNBRIGHT ME-020EA | 2,000 |
| SUNBRIGHT ME-050EA | 5,000 |
| SUNBRIGHT ME-100EA | 10,000 |
| SUNBRIGHT ME-200EA | 20,000 |
| SUNBRIGHT ME-300EA | 30,000 |
| SUNBRIGHT ME-400EA | 40,000 |

New


•SUNBRIGHT® SH Series (Thiol PEG)

| Product name | MW |
|--------------------|--------|
| SUNBRIGHT ME-020SH | 2,000 |
| SUNBRIGHT ME-050SH | 5,000 |
| SUNBRIGHT ME-100SH | 10,000 |
| SUNBRIGHT ME-200SH | 20,000 |
| SUNBRIGHT ME-300SH | 30,000 |
| SUNBRIGHT ME-400SH | 40,000 |

•SUNBRIGHT® MA Series (Maleimide PEG)


C2 Type

| Product name | MW |
|----------------------|--------|
| SUNBRIGHT ME-020MA | 2,000 |
| SUNBRIGHT ME-050MA | 5,000 |
| SUNBRIGHT ME-100MA | 10,000 |
| SUNBRIGHT ME-200MA0B | 20,000 |
| SUNBRIGHT ME-300MA | 30,000 |
| SUNBRIGHT ME-400MA | 40,000 |

New

※ Patent No. : US6875841

C5 Type **New**

| Product name | MW |
|----------------------|--------|
| SUNBRIGHT ME-050MA3* | 5,000 |
| SUNBRIGHT ME-100MA3* | 10,000 |
| SUNBRIGHT ME-200MA3* | 20,000 |
| SUNBRIGHT ME-400MA3* | 40,000 |

New


New

New

New


* : make-to-order

•SUNBRIGHT® IA series (Iodoacetamide PEG)


| Product name | MW |
|--------------------|--------|
| SUNBRIGHT ME-200IA | 20,000 |
| SUNBRIGHT ME-300IA | 30,000 |
| SUNBRIGHT ME-400IA | 40,000 |

•SUNBRIGHT® CA series (Aminoxy PEG)

| Product name | MW |
|--------------------|--------|
| SUNBRIGHT ME-100CA | 10,000 |
| SUNBRIGHT ME-200CA | 20,000 |
| SUNBRIGHT ME-300CA | 30,000 |
| SUNBRIGHT ME-400CA | 40,000 |

•SUNBRIGHT® HZ series (Hydrazide PEG)

| Product name | MW |
|---------------------|--------|
| SUNBRIGHT ME-200HZ* | 20,000 |
| SUNBRIGHT ME-300HZ* | 30,000 |
| SUNBRIGHT ME-400HZ* | 40,000 |

* : make-to-order

2. Bi-Functional PEGs

Bi-functional PEGs are the most popular derivatives for crosslinking proteins, enzymes and other pharmaceutical substance. Recently, these PEGs are applied for hydrogel components for base matrix of iPS cells or EPS cells in tissue engineering.

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| Product name | X | MW |
|---------------------|--|--------|
| NHS PEG | | |
| SUNBRIGHT DE-034CS* | Succinimidylsuccinate -CO-CH ₂ CH ₂ -COO-NHS | 3,400 |
| SUNBRIGHT DE-034GS* | Succinimidylglutarate -CO-CH ₂ CH ₂ CH ₂ COO-NHS | 3,400 |
| SUNBRIGHT DE-100GS | | 10,000 |
| SUNBRIGHT DE-200GS | | 20,000 |
| SUNBRIGHT DE-030AS* | Succinimidylcarboxymethyl -CH ₂ -COO-NHS | 3,000 |
| SUNBRIGHT DE-034HS* | Succinimidylcarboxypentyl -CH ₂ CH ₂ CH ₂ CH ₂ CH ₂ -COO-NHS | 3,400 |
| SUNBRIGHT DE-100HS | | 10,000 |
| SUNBRIGHT DE-200HS | | 20,000 |
| Amine PEG | | |
| SUNBRIGHT DE-010PA | Propylamine -CH ₂ CH ₂ CH ₂ NH ₂ | 1,000 |
| SUNBRIGHT DE-020PA | | 2,000 |
| SUNBRIGHT DE-034PA* | | 3,400 |
| SUNBRIGHT DE-100PA | | 10,000 |
| SUNBRIGHT DE-200PA | | 20,000 |
| SUNBRIGHT DE-300PA | | 30,000 |
| Thiol PEG | | |
| SUNBRIGHT DE-034SH | Ethanethiol -CH ₂ CH ₂ SH | 3,400 |
| SUNBRIGHT DE-100SH* | | 10,000 |
| SUNBRIGHT DE-200SH* | | 20,000 |
| Maleimide PEG | | |
| SUNBRIGHT DE-100MA | Maleimide | 10,000 |
| SUNBRIGHT DE-200MA | -(CH ₂) ₃ -NHCO-CH ₂ CH ₂ -Maleimide | 20,000 |

* : make-to-order


3. Multi-Arm PEGs

We supply various derivatives in which functional groups are attached to the terminals of multi-arm (4-arm and 8-arm) PEGs. Multi-arm PEGs can be used not only as the starting materials for hydrogels, but also as modifiers to improve the solubility of hydrophobic drugs.

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3-1. 4-arm-Functional PEG Series

PTE Series (Pentaerythritol, tetra-polyethylene glycol ether)


| Product name | X | MW |
|------------------------------------|--|--------|
| NHS PEG | | |
| SUNBRIGHT PTE-100CS | Succinimidylsuccinate -CO-CH ₂ CH ₂ -COO-NHS | 10,000 |
| SUNBRIGHT PTE-050GS* | | 5,000 |
| SUNBRIGHT PTE-100GS | | 10,000 |
| SUNBRIGHT PTE-150GS* | Succinimidylglutarate -CO-CH ₂ CH ₂ CH ₂ -COO-NHS | 15,000 |
| SUNBRIGHT PTE-200GS | | 20,000 |
| SUNBRIGHT PTE-400GS* | | 40,000 |
| SUNBRIGHT PTE-100HS* | | 10,000 |
| SUNBRIGHT PTE-200HS* | Succinimidylcarboxypentyl -CH ₂ CH ₂ CH ₂ CH ₂ CH ₂ -COO-NHS | 20,000 |
| SUNBRIGHT PTE-400HS* | | 40,000 |
| p-Nitrophenyl carbonate PEG | | |
| SUNBRIGHT PTE-100NP | | 10,000 |
| SUNBRIGHT PTE-200NP | p-Nitrophenylcarbonate -COO-phenyl-NO ₂ | 20,000 |
| SUNBRIGHT PTE-400NP* | | 40,000 |
| Amine PEG | | |
| SUNBRIGHT PTE-100PA | | 10,000 |
| SUNBRIGHT PTE-150PA* | Propylamine -CH ₂ CH ₂ CH ₂ NH ₂ | 15,000 |
| SUNBRIGHT PTE-200PA | | 20,000 |
| SUNBRIGHT PTE-400PA* | | 40,000 |
| Thiol PEG | | |
| SUNBRIGHT PTE-050SH* | Ethanethiol -CH ₂ CH ₂ SH | 5,000 |
| SUNBRIGHT PTE-100SH | | 10,000 |
| SUNBRIGHT PTE-200SH | | 20,000 |
| Maleimide PEG | | |
| SUNBRIGHT PTE-100MA* | Maleimide | 10,000 |
| SUNBRIGHT PTE-200MA* | -(CH ₂) ₃ -NHCO-CH ₂ CH ₂ -Maleimide | 20,000 |
| SUNBRIGHT PTE-400MA* | | 40,000 |

* : make-to-order

3-2. 8-arm-Functional PEG Series

HGEO Series (Hexa-glycerine, octa-polyethylene glycol ether)

| Product name | X | MW |
|------------------------------------|---|--------|
| NHS PEG | | |
| SUNBRIGHT HGEO-150CS* | Succinimidylsuccinate -CO-CH ₂ CH ₂ -COO-NHS | 15,000 |
| SUNBRIGHT HGEO-150GS* | | 15,000 |
| SUNBRIGHT HGEO-200GS* | Succinimidylglutarate -CO-CH ₂ CH ₂ CH ₂ -COO-NHS | 20,000 |
| SUNBRIGHT HGEO-400GS* | | 40,000 |
| p-Nitrophenyl carbonate PEG | | |
| SUNBRIGHT HGEO-200NP* | p-Nitrophenylcarbonate | 20,000 |
| SUNBRIGHT HGEO-400NP* | -COO-phenyl-NO ₂ | 40,000 |
| Amine PEG | | |
| SUNBRIGHT HGEO-150PA* | Propylamine -CH ₂ CH ₂ CH ₂ NH ₂ | 15,000 |
| Thiol PEG | | |
| SUNBRIGHT HGEO-200SH* | Ethanethiol -CH ₂ CH ₂ SH | 20,000 |

* : make-to-order

4. Branched PEGs



NOF provides 2-branched Activated PEGs, as illustrated below. Both the molecular weights of the respective PEG chains and the linker moieties leading to the functional groups are adjustable. These are available with amino, aldehyde, maleimide, carbonate and carboxylate activated esters as the terminal functional groups.

Furthermore, NOF has also developed newer types of Branched Activated PEGs, such as 3-branched and 4-branched derivatives. Although aqueous solutions of higher-molecular PEG derivatives tend to show increased viscosity, use of higher-branched Activated PEGs, such as 3-branched and 4-branched derivatives, reduces the viscosity of the relevant aqueous solutions, even without altering the molecular weights.

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
4-1. 2-arm Branched PEG (SUNBRIGHT® GL2 series)

Total Mw 40kDa: 20k × 2
60kDa: 30k × 2
80kDa: 40k × 2


GS2 Type (NHS-ester PEG)

| Product name | MW |
|-----------------------|--------|
| SUNBRIGHT GL2-200GS2 | 20,000 |
| SUNBRIGHT GL2-400GS2 | 40,000 |
| SUNBRIGHT GL2-600GS2 | 60,000 |
| SUNBRIGHT GL2-800GS2* | 80,000 |

* : make-to-order


TS Type (NHS Carbonate PEG)

| Product name | MW |
|----------------------|--------|
| SUNBRIGHT GL2-200TS | 20,000 |
| SUNBRIGHT GL2-400TS | 40,000 |
| SUNBRIGHT GL2-600TS | 60,000 |
| SUNBRIGHT GL2-800TS* | 80,000 |

* : make-to-order


NP Type (P-Nitrophenyl Carbonate PEG)

| Product name | MW |
|----------------------|--------|
| SUNBRIGHT GL2-100NP | 10,000 |
| SUNBRIGHT GL2-200NP | 20,000 |
| SUNBRIGHT GL2-400NP | 40,000 |
| SUNBRIGHT GL2-600NP | 60,000 |
| SUNBRIGHT GL2-800NP* | 80,000 |

* : make-to-order


AL3 Type (Aldehyde PEG)

| Product name | MW |
|-----------------------|--------|
| SUNBRIGHT GL2-200AL3 | 20,000 |
| SUNBRIGHT GL2-400AL3 | 40,000 |
| SUNBRIGHT GL2-600AL3 | 60,000 |
| SUNBRIGHT GL2-800AL3* | 80,000 |


*: make-to-order

PA Type (Aminopropyl PEG)


| Product name | MW |
|---------------------|--------|
| SUNBRIGHT GL2-200PA | 20,000 |
| SUNBRIGHT GL2-400PA | 40,000 |
| SUNBRIGHT GL2-600PA | 60,000 |
| SUNBRIGHT GL2-800PA | 80,000 |

C2 MA Type (Maleimide PEG)

| Product name | MW |
|---------------------|--------|
| SUNBRIGHT GL2-200MA | 20,000 |
| SUNBRIGHT GL2-400MA | 40,000 |
| SUNBRIGHT GL2-600MA | 60,000 |
| SUNBRIGHT GL2-800MA | 80,000 |


C5 MA Type (Maleimide PEG) **New**

| Product name | MW |
|-----------------------|--------|
| SUNBRIGHT GL2-400MA3* | 40,000 |

New

*: make-to-order



CA Type (Aminoxy PEG)

| Product name | MW |
|----------------------|--------|
| SUNBRIGHT GL2-100CA* | 10,000 |
| SUNBRIGHT GL2-200CA | 20,000 |
| SUNBRIGHT GL2-400CA | 40,000 |
| SUNBRIGHT GL2-600CA* | 60,000 |
| SUNBRIGHT GL2-800CA* | 80,000 |

*: make-to-order

4-2. 3-arm Branched PEG (SUNBRIGHT® GL3 series)



Total Mw 50kDa
Mw of (CH₂CH₂O)_m part :about 10kDa

| Product name | X | MW |
|------------------------------------|--|--------|
| NHS PEG | | |
| SUNBRIGHT GL3-400GS100U | Succinimidylglutarate -CO-CH ₂ CH ₂ CH ₂ -COO-NHS | 50,000 |
| SUNBRIGHT GL3-400HS100U | Succinimidylcarboxypentyl -CH ₂ CH ₂ CH ₂ CH ₂ CH ₂ -COO-NHS | 50,000 |
| SUNBRIGHT GL3-400TS100U | Succinimidylcarboxy -COO-NHS | 50,000 |
| p-Nitrophenyl carbonate PEG | | |
| SUNBRIGHT GL3-400NP100U | p-Nitrophenylcarbonate -COO-phenyl-NO ₂ | 50,000 |
| Aldehyde PEG | | |
| SUNBRIGHT GL3-400AL100U | Aldehyde -CH ₂ -CH ₂ -CHO | 50,000 |
| Amine PEG | | |
| SUNBRIGHT GL3-400PA100U | Propylamine -CH ₂ CH ₂ CH ₂ NH ₂ | 50,000 |
| Maleimide PEG | | |
| SUNBRIGHT GL3-400MA100U | Maleimide -(CH ₂) ₃ -NHCO-CH ₂ CH ₂ -maleimide | 50,000 |
| Aminoxy PEG | | |
| SUNBRIGHT GL3-400CA100U | Aminoxy -CONH-CH ₂ CH ₂ -ONH ₂ | 50,000 |

4-3. 4-arm Branched PEG

1) SUNBRIGHT® GL4 series


Total Mw 40kDa: $5\text{k} \times 2, 7.5\text{k} \times 4$
 60kDa: $5\text{k} \times 2, 12.5\text{k} \times 4$
 80kDa: $5\text{k} \times 2, 17.5\text{k} \times 4$

| Product name | X | MW |
|-----------------------------|---|--------|
| NHS PEG | | |
| SUNBRIGHT GL4-400GS2 | Succinimidylglutarate-aminopropoxy -CH ₂ CH ₂ CH ₂ -NHCO-(CH ₂) ₃ -COO-NHS | 40,000 |
| SUNBRIGHT GL4-600GS2* | | 60,000 |
| SUNBRIGHT GL4-800GS2* | | 80,000 |
| SUNBRIGHT GL4-400TS | Succinimidylcarboxy -COO-NHS | 40,000 |
| SUNBRIGHT GL4-600TS* | | 60,000 |
| SUNBRIGHT GL4-800TS* | | 80,000 |
| p-Nitrophenyl carbonate PEG | | |
| SUNBRIGHT GL4-400NP | p-Nitrophenylcarbonate -COO-phenyl-NO ₂ | 40,000 |
| SUNBRIGHT GL4-600NP* | | 60,000 |
| SUNBRIGHT GL4-800NP* | | 80,000 |
| Aldehyde PEG | | |
| SUNBRIGHT GL4-400AL3 | Aldehyde -CONH-CH ₂ CH ₂ -CHO | 40,000 |
| SUNBRIGHT GL4-600AL3* | | 60,000 |
| SUNBRIGHT GL4-800AL3* | | 80,000 |
| Amine PEG | | |
| SUNBRIGHT GL4-400PA | Propylamine -CH ₂ CH ₂ CH ₂ NH ₂ | 40,000 |
| SUNBRIGHT GL4-600PA* | | 60,000 |
| SUNBRIGHT GL4-800PA* | | 80,000 |
| Maleimide PEG | | |
| SUNBRIGHT GL4-400MA | Maleimide -(CH ₂) ₃ -NHCO-CH ₂ CH ₂ -maleimide | 40,000 |
| SUNBRIGHT GL4-600MA* | | 60,000 |
| SUNBRIGHT GL4-800MA* | | 80,000 |
| Aminoxy PEG | | |
| SUNBRIGHT GL4-400CA | Aminoxy -CONH-CH ₂ CH ₂ -ONH ₂ | 40,000 |
| SUNBRIGHT GL4-600CA* | | 60,000 |
| SUNBRIGHT GL4-800CA* | | 80,000 |

* : make-to-order

2) SUNBRIGHT® XY4 series New

| Product name | X | MW | |
|------------------------------------|---|--------|--------------------------------------|
| NHS PEG | | | |
| SUNBRIGHT XY4-200GS2* | Succinimidylglutarate-aminopropoxy | 20,000 | New |
| SUNBRIGHT XY4-400GS2 | -CH ₂ CH ₂ CH ₂ -NHCO-(CH ₂) ₃ -COO-NHS | 40,000 | New |
| SUNBRIGHT XY4-200TS* | Succinimidylcarboxy | 20,000 | New |
| SUNBRIGHT XY4-400TS | -COO-NHS | 40,000 | New |
| p-Nitrophenyl carbonate PEG | | | |
| SUNBRIGHT XY4-200NP* | p-Nitrophenylcarbonate | 20,000 | New |
| SUNBRIGHT XY4-400NP | -COO-phenyl-NO ₂ | 40,000 | New |
| Aldehyde PEG | | | |
| SUNBRIGHT XY4-200AL3* | Aldehyde | 20,000 | New |
| SUNBRIGHT XY4-400AL3 | -CONH-CH ₂ CH ₂ -CHO | 40,000 | New |
| Maleimide PEG | | | |
| SUNBRIGHT XY4-200MA* | Maleimide | 20,000 | New |
| SUNBRIGHT XY4-400MA | -(CH ₂) ₃ -NHCO-CH ₂ CH ₂ -maleimide | 40,000 | New |

* : make-to-order

Comparison Data of PEG-Solution Viscosity


| Compound | Structure | Viscosity (mPa · s) |
|-----------------|-----------|---------------------|
| MeO-PEG (40kDa) | Linear | 19.3 |
| GL2 (40kDa) | 2-arm | 18.8 |
| GL4 (40kDa) | 4-arm | 10.9 |
| XY4 (40kDa) | 4-arm | 9.6 |

※ as aqueous solution at 10 wt%, at 40°C

4-4. Lysine Branched PEG New (SUNBRIGHT® LY series)

Collaborating with Enzon Pharmaceuticals, Inc., NOF has a right to manufacture and sell their unique Lysine Branched PEG (NOF's Trade Name is SUNBRIGHT LY series) for customer research purpose only. Lysine Branched PEG has been already used in several marketed PEGylated-drugs such as PEGASIS®, MACUGEN® and CIMZIA® for 10 years.

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* : make-to-order

※ The Lysine Branched PEG belong to ENZON Pharmaceuticals, Inc., is sold by NOF for research and development purpose only.


5. Heterofunctional PEGs

X-(OCH₂CH₂)_n-Y

By using hetero-type activated PEGs, different molecules can be conjugated onto the each end of the PEGs. They are also useful for surface modification.

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
•Hydroxy-PEG-Amine

| Product name | MW |
|---------------------|--------|
| SUNBRIGHT HO-020PA | 2,000 |
| SUNBRIGHT HO-034PA | 3,400 |
| SUNBRIGHT HO-050PA | 5,000 |
| SUNBRIGHT HO-100PA | 10,000 |
| SUNBRIGHT HO-200PA | 20,000 |
| SUNBRIGHT HO-400PA* | 40,000 |

* : make-to-order


•Amino-PEG-Carboxylic acid

| Product name | MW |
|---------------------|-------|
| SUNBRIGHT PA-020HC* | 2,000 |
| SUNBRIGHT PA-034HC* | 3,400 |
| SUNBRIGHT PA-050HC* | 5,000 |

* : make-to-order

•Boc-protected-Amino-PEG-Carbonate-NHS


| Product name | MW |
|---------------------|--------|
| SUNBRIGHT BO-020TS* | 2,000 |
| SUNBRIGHT BO-034TS* | 3,400 |
| SUNBRIGHT BO-050TS* | 5,000 |
| SUNBRIGHT BO-100TS* | 10,000 |
| SUNBRIGHT BO-200TS* | 20,000 |

New

* : make-to-order

•Maleimide-PEG-Carbonate-NHS


| Product name | MW |
|---------------------|--------|
| SUNBRIGHT MA-020TS* | 2,000 |
| SUNBRIGHT MA-034TS* | 3,400 |
| SUNBRIGHT MA-050TS* | 5,000 |
| SUNBRIGHT MA-100TS* | 10,000 |
| SUNBRIGHT MA-200TS* | 20,000 |
| SUNBRIGHT MA-400TS* | 40,000 |

New


* : make-to-order

•Hydroxy-PEG-Aldehyde

| Product name | MW |
|--------------------|--------|
| SUNBRIGHT HO-050AL | 5,000 |
| SUNBRIGHT HO-100AL | 10,000 |
| SUNBRIGHT HO-200AL | 20,000 |
| SUNBRIGHT HO-300AL | 30,000 |

•Biotin-PEG-Carbonate-NHS

| Product name | MW |
|---------------------|-------|
| SUNBRIGHT BI-050TS* | 5,000 |

* : make-to-order

6. Forked PEGs New


•SUNBRIGHT® PTE2 Series

Forked structures have the advantage of placing two reactive groups at precise distances apart. These "Forked" PEGs have become very popular for mimicking the heavy chain domain in an antibody or fragment antibody and other applications where two proteins held in proximity are advantageous.

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7. Releasable PEGs New

•SUNBRIGHT® BE-Linker series


Pegylated therapeutic protein or small molecule pharmacology would be substantially expanded if the original therapeutic proteins or small molecules could be regenerated in vivo. Diminution of activity of both enzymes and protein ligands is commonly encountered following permanent conjugation with polyethylene glycol (PEG). However, Enzon's releasable PEG linkers have the ability to improve the pharmaceutical properties of therapeutic proteins or small molecules through releasable PEGylation and maintain the therapeutic effect of these molecules by regenerating the bioactive proteins or small molecules in vivo. To make available Enzon's releasable linker technology, NOF introduces the Benzyl Elimination (BE) linkers, i.e. PEG-BE-NHS, SUNBRIGHT BE Series, as listed in the following.

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The advantage of SUNBRIGHT® BE series:

- Controls pharmacokinetics of drugs in vivo (C_{max} , T_{max} , $T_{1/2}$)
- Provides high loading capability
- Deposits poorly soluble drugs at target tissue
- Customized cellular uptake, trafficking and processing

Releasing drug from PEG-drug conjugates (NOF SUNBRIGHT® BE-Linker series)

| Product name | Chemical structure | MW | |
|---------------------------|--------------------|--------|--------------------------------------|
| SUNBRIGHT ME-BE200CH-TS* | | 20,000 | New |
| SUNBRIGHT ME-BE200CM-TS* | | 20,000 | New |
| SUNBRIGHT ME-BE200EH5-TS* | | 20,000 | New |
| SUNBRIGHT ME-BE200EM5-TS* | | 20,000 | New |


※ The Releasable Activated PEG belong to ENZON Pharmaceuticals Inc., is sold by NOF for research and development purpose only.

※ Orders for other molecular weights and linkers may require custom synthesis. Please e-mail at ddsinfo@nof.co.jp for a quote on custom reagents.

* : make-to-order


Characteristics of SUNBRIGHT® BE-Linker

•Effects of PEG-BE Linker on Lysozyme Regeneration

Regeneration of lysozyme from PEG-BE-lysozyme conjugates in rat plasma.
S Lee, et al., *Bioconjugate Chem.* **12(2)**, 163-169 (2001)

•Antitumor effects of PEG-BE-immunotoxin (BE-PEG-SS1P)

Antitumor effects of PEG-BE-SS1P and SS1P on A431-K5 solid tumors.
Cells (3×10^6) were injected s.c. into nude mice on day 0. On day 7 (tumor volume = $\sim 140 \text{ mm}^3$), animals were treated with i.v. injections of the compounds at the indicated doses. Control mice received vehicle only.
D Filpula, et al., *Bioconjugate Chem.* **18(3)**, 773-784 (2007)

•Hydrolysis, Cytotoxicity, Antitumor activity of PEG-BE-Daunorubicin

| Compound | $t_{1/2}(\text{hr})$ buffer (pH 7.4) | $t_{1/2}(\text{hr})$ rat plasma | IC ₅₀ * ¹ (nM) p388/0 | Antitumor * ² i.v. % T/C * ³ |
|--|--|---------------------------------------|---|--|
| Daunorubicin (DNR) | control | control | 3 | 117.0 |
| Esters $\text{PEG}-\text{C}(=\text{O})-\text{O}-\text{C}_6\text{H}_3(\text{R}_1, \text{R}_2)-\text{CH}_2\text{O}-\text{C}(=\text{O})-\text{NH}-\text{DNR}$ $\text{R}_1 = \text{R}_2 = \text{H}$ $\text{R}_1 = \text{R}_2 = \text{OCH}_3$ $\text{R}_1 = \text{R}_2 = \text{CH}_3$ | >24 | 0.4 | 8 | NA |
| | >48 | 1.0 | 27 | 48.2 |
| | >48 | 1.9 | 55 | 67.9 |
| Carbonates $\text{PEG}-\text{O}-\text{C}(=\text{O})-\text{O}-\text{C}_6\text{H}_3(\text{CH}_3)_2-\text{CH}_2\text{O}-\text{C}(=\text{O})-\text{NH}-\text{DNR}$ | >48 | 2.9 | 179 | 74.4 |

*¹ : Half maximal (50%) inhibitory concentration, Inhibition of P388 leukemia

*² : 3 mg/kg/dose of DNR to BALB/c mice bearing subcutaneous Madison lung carcinoma(M-109) on days 1 and 4(ip) / 3 and 6(iv).


*³ : Percent treatment over control (%T/C)

RB Greenwald et al., *J. Med. Chem.* **42**, 3657-3667 (1999)

8. Comb-shaped co-polymers (carboxylic anhydride type)

When these reactive polymers are used for modification of enzymes, less reactivity is required to introduce more PEG chains onto the enzyme surfaces than that of a single polymeric PEG modifier. Consequently, the activities of the modified enzymes can be improved.

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
General Characteristics of PEG SUNBRIGHT® AM series

| Product Name | MW of AO (PEG) | DP ⁽¹⁾ | Total MW |
|--------------------|----------------|-------------------|---------------------|
| SUNBRIGHT AM-0530K | ca. 500 | ca. 30 ~ 40 | ca. 15,000 ~ 20,000 |
| SUNBRIGHT AM-1510K | ca. 1,500 | ca. 10 ~ 15 | ca. 15,000 ~ 20,000 |
| SUNBRIGHT AM-2090P | ca. 2,000 | ca. 10 ~ 20 | ca. 20,000 ~ 40,000 |


(1) Degree of Polymerization


Moreover, with the use of this compound, the hydrophilic/lipophilic balance can be adjusted by changing the composition ratios of the ethylene oxide and propylene oxide in the polyalkylene glycol present in this compound.

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Efficacy of Comb-type PEG Enzyme Conjugation Stability of Modified Protease in water (at 40°C)

Stability of modified L-Asparaginase

References: for modified enzymes using SUNBRIGHT AM-series:

- 1) J. Kajiuchi et al., *J. Chem. Eng. Japan* **25**, 202 (1992)
- 2) T. Masunaga, et al., *J. SCCJ* **27(3)**, 276 (1993)
- 3) H. Sasaki et al., *Biochem. Biophys. Res. Commun* **197**, 287 (1993)
- 4) Y. Kodera, et al., *Bioconjugate Chem.* **5**, 283 (1994)
- 5) K. Matsuo, et al., *Adv. Bioseparation Eng.* 1993 **3**, 56 (1994)

III

Phospholipids and Lipids for Liposomal Formulations <COATSOME® Series>

Using the latest synthetic and purification technologies, NOF manufactures and supplies highly purified phospholipid derivatives (COATSOME Series) suitable for lipid emulsion formulations and liposome formulations which are extensively used for pharmaceuticals and cosmetics. Various kinds and grades of such phospholipids are readily available. In addition, we can respond to requests from our customers with tailor-made synthesis of any phospholipid with the desired structures. To benefit easy liposomal capture of drugs at customer facilities, NOF provides empty liposome kits, liposomal formulations and liposomal drug delivery systems.


When PEGylated phospholipids are used as liposomes, the aqueous corona at the liposome surface facilitates stable dispersion in aqueous solutions.

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Liposomes were imaged by FEI Ltd., using Cryo-TEM (TECNAI F20TWIN Cryo).


NOF Phospholipid Product Summary

1. Purified Phospholipids from natural sources

1-1. Soybean Phospholipids

NOF manufactures high-purity lecithin containing higher contents of phosphatidylcholine by using strictly selected natural soybeans uncontaminated with genetically modified soybeans. The high-purity product is useful to obtain favorable color and odor of the final product. What is more, the less unsaturated double-bonds in the acyl chains from higher hydrogenation result in excellent stability, together with easier handling of this lecithin in its powder form.

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
| | Product Name | Description | PC purity | MW |
|-------------|-----------------|--|-----------|--------|
| HSPC | COATSOME NC-21E | Hydrogenated Soybean phosphatidylcholine | 98% up | Ca.785 |
| HSPC | COATSOME NC-21 | Hydrogenated Soybean phosphatidylcholine | 90% up | Ca.785 |

1-2. Egg yolk Phospholipids

Egg yolk phospholipids are manufactured using virusfree yolk as the starting material, followed by removal of impurities through strict purification procedures. Unique NOF advantages include major phospholipids containing two different acyl groups (saturated and unsaturated fatty acids) in the same molecule. These unique chemical structures provide excellent dispersing ability in aqueous solution and favorable surfactant effects. Egg yolk phospholipids are also employed as carriers for lipophilic drugs as one of the major ingredients of lipid microspheres.

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
| | Product Name | Description | PC purity | MW |
|-------------|----------------|--|-----------|--------|
| EPC | COATSOME NC-50 | Non-hydrogenated Egg phosphatidylcholine | 95% up | Ca.773 |
| HEPC | COATSOME NC-11 | Hydrogenated Egg phosphatidylcholine | 95% up | Ca.777 |

1-3. Sphingomyelin

Phospholipids are known to be the primary constituents of biological membranes. Recently, however, it has been clarified that phospholipids are not only constituents of the cell membrane, but also play important direct roles in the cell membrane functions. Among the phospholipids, sphingomyelin has been revealed to be significantly important in both the formation and maintenance of lipid rafts. Considering that the lipid rafts have been shown to be involved in various signaling pathways, including immunological responses and transportation of specific materials, sphingomyelin has drawn much attention as a substance with important roles in the expression of specific cellular functions, such as intracellular information transmission and maintenance of membrane structure. Despite containing long fatty acid chains, sphingomyelin exhibits some advantageous characteristics such as a low phase-transition temperatures (T_c), that render it uniquely suitable for liposomal formations. NOF supplies high-purity sphingomyelin derived from milk and egg yolk. Concerning yolk-derived sphingomyelin, palmitic acid accounts for about 80% of the fatty acid chains bound by amide bonds; in sharp contrast, however, milk-derived sphingomyelin contains a broader range of long-chain fatty acids, including palmitic acid or longer-chain fatty acids. Because of these differences in the physical properties, customers can select between the two sources depending on their needs.

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
2. Highly-purified Synthetic Phospholipids

We use fatty acids that are highly purified by our own technologies for the manufacture of a variety of phospholipid derivatives; that is, NOF can synthesize derivatives containing a variety of fatty acids. Our phospholipids are manufactured in FDA-inspected facilities meeting GMP standards. Our quality-assured products can be safely used for the production of pharmaceuticals. NOF as a contracted manufacturer also synthesizes novel phospholipid derivatives and has completed DMF registrations for numerous novel derivatives; we are in a position to supply relevant lipids from test reagent grade to commercial production scales that are approved for manufacture of pharmaceuticals after clinical studies.

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
2-1. Phosphatidylcholine


| | Product name | Description | PC purity | MW | R ₁ COO | R ₂ COO | Tc |
|---------------|--------------------|---|-----------|-------|--------------------|--------------------|-----|
| DDPC | COATSOME MC-1010 | 1,2-Didecanoyl-sn-glycero-3-phosphocholine | 99% up | 565.7 | C10:0 | C10:0 | -6 |
| DLPC | COATSOME MC-2020 | 1,2-Dilauroyl-sn-glycero-3-phosphocholine | 99% up | 621.8 | C12:0 | C12:0 | 0 |
| DMPC | COATSOME MC-4040 | 1,2-Dimyristoyl-sn-glycero-3-phosphocholine | 99% up | 677.9 | C14:0 | C14:0 | 23 |
| DPPC | COATSOME MC-6060 | 1,2-Dipalmitoyl-sn-glycero-3-phosphocholine | 99% up | 734.0 | C16:0 | C16:0 | 41 |
| DSPC | COATSOME MC-8080 | 1,2-Distearoyl-sn-glycero-3-phosphocholine | 99% up | 790.2 | C18:0 | C18:0 | 55 |
| DOPC | COATSOME MC-8181 | 1,2-Dioleoyl-sn-glycero-3-phosphocholine | 99% up | 786.1 | C18:1 | C18:1 | -22 |
| DLoPC | COATSOME MC-8282 | 1,2-Dilinoleoyl-sn-glycero-3-phosphocholine | 99% up | 782.2 | C18:2 | C18:2 | -53 |
| DEPC | COATSOME MC-2121AL | 1,2-Dierucoyl-sn-glycero-3-phosphocholine | 99% up | 898.3 | C22:1 | C22:1 | -13 |
| EPA-PC | COATSOME MC-1515AL | 1,2-Dieicosapentaenoyl-sn-glycero-3-phosphocholine | 98% up | 826.2 | C20:5 | C20:5 | - |
| DHA-PC | COATSOME MC-2626AL | 1,2-Didocosahexaenoyl-sn-glycero-3-phosphocholine | 98% up | 878.2 | C22:6 | C22:6 | - |
| MPPC | COATSOME MC-4060 | 1-Myristoyl-2-palmitoyl-sn-glycero-3-phosphocholine | 99% up | 706.0 | C14:0 | C16:0 | 35 |
| MSPC | COATSOME MC-4080 | 1-Myristoyl-2-stearoyl-sn-glycero-3-phosphocholine | 99% up | 734.0 | C14:0 | C18:0 | 40 |
| PMPC | COATSOME MC-6040 | 1-Palmitoyl-2-myristoyl-sn-glycero-3-phosphocholine | 99% up | 706.0 | C16:0 | C14:0 | 28 |
| PSPC | COATSOME MC-6080 | 1-Palmitoyl-2-stearoyl-sn-glycero-3-phosphocholine | 99% up | 762.1 | C16:0 | C18:0 | 49 |
| SMPC | COATSOME MC-8040 | 1-Stearoyl-2-myristoyl-sn-glycero-3-phosphocholine | 99% up | 734.0 | C18:0 | C14:0 | 30 |
| SPPC | COATSOME MC-8060 | 1-Stearoyl-2-palmitoyl-sn-glycero-3-phosphocholine | 99% up | 762.1 | C18:0 | C16:0 | 44 |
| MOPC | COATSOME MC-4081 | 1-Myristoyl-2-oleoyl-sn-glycero-3-phosphocholine | 99% up | 732.0 | C14:0 | C18:1 | - |
| POPC | COATSOME MC-6081 | 1-Palmitoyl-2-oleoyl-sn-glycero-3-phosphocholine | 99% up | 760.1 | C16:0 | C18:1 | -3 |
| SOPC | COATSOME MC-8081 | 1-Stearoyl-2-oleoyl-sn-glycero-3-phosphocholine | 99% up | 788.1 | C18:0 | C18:1 | 6 |

2-2. Lyso-Phosphatidylcholine


| | Product name | Description | PC purity | MW | R ₁ COO |
|-----------------|-----------------|--|-----------|-------|--------------------|
| M-LysoPC | COATSOME MC-40H | 1-Myristoyl-2-lyso-sn-glycero-3-phosphocholine | 99% up | 467.6 | C14:0 |
| P-LysoPC | COATSOME MC-60H | 1-Palmitoyl-2-lyso-sn-glycero-3-phosphocholine | 99% up | 495.6 | C16:0 |
| S-LysoPC | COATSOME MC-80H | 1-Stearoyl-2-lyso-sn-glycero-3-phosphocholine | 99% up | 523.7 | C18:0 |
| O-LysoPC | COATSOME MC-81H | 1-Oleoyl-2-lyso-sn-glycero-3-phosphocholine | 99% up | 521.7 | C18:1 |

2-3. Phosphatidylglycerol


| | Product name | Description | PG purity | MW | R ₁ COO | R ₂ COO |
|----------------------------|--------------------|--|-----------|-------|--------------------|--------------------|
| HSPG-Na | COATSOME NG-21LS | Hydrogenated Soybean phosphatidylglycerol, sodium salt | 95% up | | | |
| EPG-Na | COATSOME NG-50LS | Non-hydrogenated Egg phosphatidylglycerol, sodium salt | 95% up | | | |
| DLPG-Na | COATSOME MG-2020LS | 1,2-Dilauroyl-sn-glycero-3-phosphoglycerol, sodium salt | 99% up | 632.8 | C12:0 | C12:0 |
| DMPG-Na | COATSOME MG-4040LS | 1,2-Dimyristoyl-sn-glycero-3-phosphoglycerol, sodium salt | 99% up | 688.9 | C14:0 | C14:0 |
| DMPG-NH₄ | COATSOME MG-4040LA | 1,2-Dimyristoyl-sn-glycero-3-phosphoglycerol, ammonium salt | 99% up | 683.9 | C14:0 | C14:0 |
| DPPG-Na | COATSOME MG-6060LS | 1,2-Dipalmitoyl-sn-glycero-3-phosphoglycerol, sodium salt | 99% up | 745.0 | C16:0 | C16:0 |
| DPPG-NH₄ | COATSOME MG-6060LA | 1,2-Dipalmitoyl-sn-glycero-3-phosphoglycerol, ammonium salt | 98% up | 740.0 | C16:0 | C16:0 |
| DSPG-Na | COATSOME MG-8080LS | 1,2-Distearoyl-sn-glycero-3-phosphoglycerol, sodium salt | 99% up | 801.1 | C18:0 | C18:0 |
| DSPG-NH₄ | COATSOME MG-8080LA | 1,2-Distearoyl-sn-glycero-3-phosphoglycerol, ammonium salt | 99% up | 796.1 | C18:0 | C18:0 |
| DOPG-Na | COATSOME MG-8181LS | 1,2-Dioleoyl-sn-glycero-3-phosphoglycerol, sodium salt | 98% up | 797.0 | C18:1 | C18:1 |
| DEPG-Na | COATSOME MG-2121LS | 1,2-Dierucoyl-sn-glycero-3-phosphoglycerol, sodium salt | 98% up | 909.2 | C22:1 | C22:1 |
| POPG-Na | COATSOME MG-6081LS | 1-Palmitoyl-2-oleoyl-sn-glycero-3-phosphoglycerol, sodium salt | 98% up | 771.0 | C16:0 | C18:1 |

2-4. Phosphatidic acid

| | Product name | Description | PA purity | MW | R ₁ COO | R ₂ COO |
|----------------|--------------------|---|-----------|-------|--------------------|--------------------|
| DLPA-Na | COATSOME MA-2020LS | 1,2-Dilauroyl-sn-glycero-3-phosphatidic acid, sodium salt | 98% up | 558.7 | C12:0 | C12:0 |
| DMPA-Na | COATSOME MA-4040LS | 1,2-Dimyristoyl-sn-glycero-3-phosphatidic acid, sodium salt | 98% up | 614.8 | C14:0 | C14:0 |
| DPPA-Na | COATSOME MA-6060LS | 1,2-Dipalmitoyl-sn-glycero-3-phosphatidic acid, sodium salt | 99% up | 670.9 | C16:0 | C16:0 |
| DSPA-Na | COATSOME MA-8080LS | 1,2-Distearoyl-sn-glycero-3-phosphatidic acid, sodium salt | 99% up | 727.0 | C18:0 | C18:0 |


2-5. Phosphatidylethanolamine

| | Product name | Description | PE purity | MW | R ₁ COO | R ₂ COO |
|--------------|--------------------|---|-----------|-------|--------------------|--------------------|
| HSPE | COATSOME NE-21* | Hydrogenated Soybean phosphatidylethanolamine | 95% up | | | |
| EPE | COATSOME NE-50* | Non-hydrogenated Egg phosphatidylethanolamine | 95% up | | | |
| DLPE | COATSOME ME-2020 | 1,2-Dilauroyl-sn-glycero-3-phosphoethanolamine | 99% up | 579.8 | C12:0 | C12:0 |
| DMPE | COATSOME ME-4040 | 1,2-Dimyristoyl-sn-glycero-3-phosphoethanolamine | 99% up | 635.9 | C14:0 | C14:0 |
| DPPE | COATSOME ME-6060 | 1,2-Dipalmitoyl-sn-glycero-3-phosphoethanolamine | 99% up | 692.0 | C16:0 | C16:0 |
| DSPE | COATSOME ME-8080 | 1,2-Distearoyl-sn-glycero-3-phosphoethanolamine | 99% up | 748.1 | C18:0 | C18:0 |
| DOPE | COATSOME ME-8181 | 1,2-Dioleoyl-sn-glycero-3-phosphoethanolamine | 99% up | 744.0 | C18:1 | C18:1 |
| DLoPE | COATSOME ME-8282 | 1,2-Dilinoleoyl-sn-glycero-3-phosphoethanolamine | 99% up | 740.0 | C18:2 | C18:2 |
| DEPE | COATSOME ME-2121AL | 1,2-Dierucoyl-sn-glycero-3-phosphoethanolamine | 99% up | 856.3 | C22:1 | C22:1 |
| POPE | COATSOME ME-6081 | 1-Palmitoyl-2-oleoyl-sn-glycero-3-phosphoethanolamine | 99% up | 718.0 | C16:0 | C18:1 |

* : make-to-order

2-6. Phosphatidylserine

| | Product name | Description | PS purity | MW | R ₁ COO | R ₂ COO |
|----------------|--------------------|--|-----------|-------|--------------------|--------------------|
| DMPS-Na | COATSOME MS-4040LS | 1,2-Dimyristoyl-sn-glycero-3-phospho-L-serine, sodium salt | 85% up | 701.8 | C14:0 | C14:0 |
| DPPS-Na | COATSOME MS-6060LS | 1,2-Dipalmitoyl-sn-glycero-3-phospho-L-serine, sodium salt | 85% up | 758.0 | C16:0 | C16:0 |
| DSPS-Na | COATSOME MS-8080LS | 1,2-Distearoyl-sn-glycero-3-phospho-L-serine, sodium salt | 85% up | 814.1 | C18:0 | C18:0 |
| DOPS-Na | COATSOME MS-8181LS | 1,2-Dioleoyl-sn-glycero-3-phospho-L-serine, sodium salt | 97% up | 810.0 | C18:1 | C18:1 |
| POPS-Na | COATSOME MS-6081LS | 1-Palmitoyl-2-oleoyl-sn-3-phospho-L-serine, sodium salt | 85% up | 784.0 | C16:0 | C18:1 |


3. PEGylated Lipids

PEGylated lipids can be applied to liposomal drug formulations to prolong the circulating plasma half-life of drugs. Coating of the liposomal surface with polyethylene glycol suppresses drug clearance in vivo by the reticuloendothelial system, thereby prolonging the half-life of the drug.

PEG-modified lipids can also be used as emulsifiers and stabilizers of microspheres in aqueous solutions.

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
3-1. PEG-phospholipids

PEG-phospholipids of the carbamate-linked type are resistant to hydrolysis, and therefore used for commercial production of PEGylated liposomes. As NOF employs its own high-quality Activated PEGs and phospholipids, high-quality PEG-phospholipids can be supplied.

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
NOF CORPORATION is the sole patent holder of these high-quality PEG-phospholipids.

SUNBRIGHT® *: make-to-order

※ Patent No. : US6679822


| Product name | Description | PEG MW | R1COO | R2COO |
|-------------------------|---|--------|-------|-------|
| SUNBRIGHT DSPE-020GL2U* | N-[Carbonyl-2',3'-Bis(methoxypolyethyleneglycol 2000)]-1,2-distearoyl-sn-glycero-3-phosphoethanolamine, sodium salt | 2,000 | C18:0 | C18:0 |
| SUNBRIGHT DSPE-050GL2U* | N-[Carbonyl-2',3'-Bis(methoxypolyethyleneglycol 5000)]-1,2-distearoyl-sn-glycero-3-phosphoethanolamine, sodium salt | 5,000 | C18:0 | C18:0 |

SUNBRIGHT® *: make-to-order

3-2. New hydrophilic phospholipid derivatives


As listed below, NOF provides new types of phospholipid derivatives bound to polyglycerin and multi-arm PEGs as hydrophilic polymers.

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SUNBRIGHT® *: make-to-order


In vivo Plasma Clearance Profile for PEG-based Liposomes and Polyglycerine-based Liposome

3-3. Diacylglycerol-PEG

Because of direct attachment of PEG chains to diacylglycerol, this product series is devoid of phosphorylethanolamine groups and other coupling chemistry. Accordingly, they are resistant to hydrolysis and are suitable for application to drugs that are vulnerable to ionic charges.

Concerning this product, NOF has officially registered its product patent in Japan.

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
SUNBRIGHT® * : make-to-order

3-4. Cholesterol-PEG derivatives

In addition to methoxypolyethyleneglycol, NOF provides various kinds of high-purity monoalkyl-substituted polyethyleneglycols. When hydrocarbon moieties comprising more than 8 carbons in the alkyl groups are used, several unique surface-active effects can be obtained. These derivatives can be employed not only for surface modification of liposomes, but also for PEGylation of proteins by the introduction of various functional groups to the PEG chain terminals. Usefulness of these products as solubilizers is also demonstrated.

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
SUNBRIGHT®

Serum clearance efficiency of new PEG-Lipid derivatives in Liposomes


References: K. Maruyama et al., *Bio. Pharm. Bull.* **19(10)**, 1347 (1996)

4. Functionalized Phospholipids

Use of functionalized phospholipids enables researchers to do protein lipidation, attaching various peptides or other biologics to the surfaces of lipid emulsions or liposomes that allows some desired target features to be obtained, and increase of the plasma half-life of the modified drugs with the introduction of PEG chains into the target drugs. NOF can supply various kinds of high-purity derivatives because we have the advantage of manufacturing our own pure, highly reactive PEG.

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
4-1. Activated phospholipids

Activated phospholipids include both maleimide groups reactive with -SH groups in biologics, and activated carboxylic esters which react with -NH₂ and -SH groups in proteins.

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

NHS-Phospholipid


| | Product name | Description | R ₁ C00 | R ₂ C00 | MW |
|----------|---------------------|--|--------------------|--------------------|-------|
| DMPE-NHS | COATSOME FE-4040SU5 | N-(Succinimidyl-oxy-glutaryl)-L-α-phosphatidylethanolamine, Dimyristoyl | C14:0 | C14:0 | 869.0 |
| DPPE-NHS | COATSOME FE-6060SU5 | N-(Succinimidyl-oxy-glutaryl)-L-α-phosphatidylethanolamine, Dipalmitoyl | C16:0 | C16:0 | 925.1 |
| DSPE-NHS | COATSOME FE-8080SU5 | N-(Succinimidyl-oxy-glutaryl)-L-α-phosphatidylethanolamine, Distearoyl | C18:0 | C18:0 | 981.2 |
| POPE-NHS | COATSOME FE-6081SU5 | N-(Succinimidyl-oxy-glutaryl)-L-α-phosphatidylethanolamine, 1-Palmitoyl-2-oleoyl | C16:0 | C18:1 | 951.2 |
| DOPE-NHS | COATSOME FE-8181SU5 | N-(Succinimidyl-oxy-glutaryl)-L-α-phosphatidylethanolamine, Dioleoyl | C18:1 | C18:1 | 977.2 |

Glu-Phospholipid


| | Product name | Description | R ₁ COO | R ₂ COO | MW |
|----------|--------------------|---|--------------------|--------------------|-------|
| DMPE-Glu | COATSOME FE-4040GL | N-Glutaryl-L-α-phosphatidylethanolamine, Dimyristoyl | C14:0 | C14:0 | 771.9 |
| DPPE-Glu | COATSOME FE-6060GL | N-Glutaryl-L-α-phosphatidylethanolamine, Dipalmitoyl | C16:0 | C16:0 | 828.1 |
| DSPE-Glu | COATSOME FE-8080GL | N-Glutaryl-L-α-phosphatidylethanolamine, Distearoyl | C18:0 | C18:0 | 884.2 |
| POPE-Glu | COATSOME FE-6081GL | N-Glutaryl-L-α-phosphatidylethanolamine, 1-Palmitoyl-2-oleoyl | C16:0 | C18:1 | 854.1 |
| DOPE-Glu | COATSOME FE-8181GL | N-Glutaryl-L-α-phosphatidylethanolamine, Dioleoyl | C18:1 | C18:1 | 880.1 |

PDP-Phospholipid **New**

| | Product name | Description | R ₁ COO | R ₂ COO | MW |
|----------|--------------------|---|--------------------|--------------------|------------------|
| DPPE-PDP | COATSOME FE-6060DT | N-[3-(2-Pyridinyldithio)-1-oxopropyl]-L-α-phosphatidylethanolamine, Dipalmitoyl | C16:0 | C16:0 | 911.2 New |


4-2. Activated PEG phospholipids

Phospholipid-PEG-NH₂

| | Product name | Description | PEG MW |
|--------------------------|-----------------------|---|--------|
| DSPE-PEG-NH ₂ | SUNBRIGHT DSPE-020PA | N-(aminopropyl polyethyleneglycol)carbamyl- distearoylphosphatidyl-ethanolamine | 2,000 |
| DSPE-PEG-NH ₂ | SUNBRIGHT DSPE-034PA* | N-(aminopropyl polyethyleneglycol)carbamyl- distearoylphosphatidyl-ethanolamine | 3,400 |
| DSPE-PEG-NH ₂ | SUNBRIGHT DSPE-050PA | N-(aminopropyl polyethyleneglycol)carbamyl- distearoylphosphatidyl-ethanolamine | 5,000 |

SUNBRIGHT® * : make-to-order


Phospholipid-PEG-MAL

| | Product name | Description | PEG MW |
|--------------|-----------------------|--|--------|
| DSPE-PEG-MAL | SUNBRIGHT DSPE-020MA | N-[(3-Maleimide-1-oxopropyl)aminopropyl polyethyleneglycol-carbamyl] distearoylphosphatidyl-ethanolamine | 2,000 |
| DSPE-PEG-MAL | SUNBRIGHT DSPE-034MA | N-[(3-Maleimide-1-oxopropyl)aminopropyl polyethyleneglycol-carbamyl] distearoylphosphatidyl-ethanolamine | 3,400 |
| DSPE-PEG-MAL | SUNBRIGHT DSPE-050MA* | N-[(3-Maleimide-1-oxopropyl)aminopropyl polyethyleneglycol-carbamyl] distearoylphosphatidyl-ethanolamine | 5,000 |

SUNBRIGHT® * : make-to-order


Phospholipid-PEG-NHS

| Product name | Description | PEG MW |
|---------------------|-----------------------|--------|
| DSPE-PEG-NHS | SUNBRIGHT DSPE-020GS | 2,000 |
| DSPE-PEG-NHS | SUNBRIGHT DSPE-034GS* | 3,400 |
| DSPE-PEG-NHS | SUNBRIGHT DSPE-050GS* | 5,000 |

SUNBRIGHT® *: make-to-order

Phospholipid-PEG-ALD New


| Product name | Description | PEG MW |
|---------------------|-----------------------|--------|
| DSPE-PEG-ALD | SUNBRIGHT DSPE-034AL* | 3,400 |


SUNBRIGHT® *: make-to-order

4-3. Fluorescent phospholipids

NBD-DPPE


Dansyl-DPPE


| Product name | Description | MW | R ₁ COO | R ₂ COO |
|--------------------|--------------------|--------|--------------------|--------------------|
| NBD-DPPE | COATSOME FE-6060NB | 956.3 | C16:0 | C16:0 |
| Dansyl-DPPE | COATSOME FE-6060DA | 1026.5 | C16:0 | C16:0 |

NBD-DPPE


Dansyl-DPPE

5. Novel Lipids and Cationic Lipids


5-1. Ether phospholipid

| Product name | Description | MW |
|------------------|------------------|--|
| C10:0 Diether PC | COATSOME EC-1010 | 1,2-Di-O-Decyl-sn-glycero-3-phosphocholine |
| | | 537.8 |

COATSOME®


5-2. DOTAP

| Product name | Description | MW |
|--------------|--------------------|--|
| DOTAP | COATSOME CL-8181TA | 1,2-Dioleoyloxy-3-trimethylammonium propane chloride |
| | | 698.5 |

COATSOME®

5-3. SAINT™-Solid-1(SN-010) New


| Product name | Description | MW |
|------------------------|--|-------|
| SAINT™-Solid-1(SN-010) | 4-[(9Z,28Z)-heptatriaconta-9,28-dien-19-yl]-1-methylpyridin-1-ium chloride | 644.5 |

New

Patent No. :EP755924, US5853694, US6726894. Licensing information is available from Synvolux Therapeutics B.V.


DNA Delivery

Comparison of GFP Expression in Murine Embryonic Stem Cells

siRNA Delivery


SAINT-RED™ Comparison
Silencing of EGP-2 in SKOV-3 Cells
[siRNA] = 40ng

Protein Delivery

Delivery of DNA Modifying Enzymes

Delivery of DNA methyltransferase M.Sss I into high E-cadherin expressing T6HBE cells

References :

- Bernardia T.F. van der Gun, et al., *J Control Release* **123**, 228-238 (2007)
 Bernardia T.F. van der Gun, et al., *Int. J. Cancer* **123**, 484-489 (2008)
 Sigridur A. Asgeisdottir, et al., *J Control Release* **141**, 241-251 (2010)
 Joanna E. Adrian, et al., *J Control Release* **144**, 341-349 (2007)

Note :

This product is covered by one or more patent applications and/or foreign counter-part patent applications owned by Synvolux Therapeutics B.V.

6. Empty Liposomes <COATSOME® EL Series>


In general, preparation of liposomal drugs requires tedious adjustments of the liposomal composition and also control of the particle size and encapsulation efficiency, which make researchers usually reluctant to formulate liposomal drugs. NOF has successfully overcome these problems by developing unique experimental kits for the investigation of liposomal drug delivery systems (DDS).

Freeze-dried ready-to-use liposome powder, called Empty Liposomes, filled in vials is composed of several kinds of phospholipids and electrolytes. When a drug solution is injected into a vial and gently shaken, the drug is easily encapsulated in the liposomes. Empty Liposomes have a special advantage in that even drugs such as the anthracyclines and aminoglycosides can be efficiently encapsulated without the need for special technologies such as extrusion.

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How to use

- (1) COATSOME EL removed from cool storage and allowed to come to room temperature.
- (2) Add aqueous solution of a drug in the temperature range of 16°C to 40°C.
- (3) Shake the vial gently three to five times by hand.
- (4) The liposomal drug is now ready for use.

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6-1. EL Series Product Characteristics


Product compositions and characteristics are shown below. Cationic charged liposomes are prepared by incorporating stearylamine as a charged lipid into the formulation. Nonionic and anionic charged liposomes are prepared by altering the ratio of DPPG as the charged lipid. While preparing an aqueous solution of the drug in water (Distilled Water for Injection), it is necessary to decide the concentrations according to the molarity of the DPPG as the charged lipid. For further details, please refer to "Suggestions for use" in the next section.

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
•COATSOME® EL-01-D

We also provide new cationic liposomes (COATSOME EL-01-D) applicable to gene delivery systems. This novel cationic liposome contains not only O, O'-ditetradecanoyl-N-(α -trimethylammonioacetyl) diethanolamine chloride (DC-6-14) as a cationic lipid, but also DOPE and cholesterol, and confers the characteristics of efficient transfection activity and expression ability both in vitro serum-containing media and in vivo assay systems.

O,O'-ditetradecanoyl- N-(α -trimethylammonioacetyl)diethanolamine chloride (DC-6-14)

| Product Name | Lipid Composition($\mu\text{mol/vial}$) | Total amount of lipids(mg) |
|------------------|---|----------------------------|
| COATSOME EL-01-D | DOPE:Cholesterol:DC-6-14=0.75 : 0.75 : 1.00 | 1.51 |

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COATSOME® EL

How to use

- When 2 mL of distilled water is added, the osmotic pressure ratio of the product becomes 0.8 to 1.1. The osmotic pressure should be adjusted according to the need in the intended experiment.
- For enhancement of the encapsulation efficiency of the drug by electrostatic interaction, the molar ratio of the charged lipids (either stearylamine or DPPG) against the drug solution plays an important role; the molar ratio should be at least more than 2, and preferably not less than 3.
- Adjust the temperature range for storage of the aqueous drug solution between 16°C and 40°C.
- Use the liposomal drugs immediately after the adjustment. (Note: Not for use in humans.)

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6-3. Determination of liposome encapsulation efficiency

Method for removal of the unencapsulated drugs

(a) Dialysis

Place the prepared drug liposome solution in a dialyzing tube and dialyze it against an isotonic solution: the unencapsulated drug leaks into the isotonic solution outside the tube.

(b) Gel filtration

Pass the prepared liposome solution through a column filled with a packing material for gel filtration. The liposomes and the unencapsulated drug are eluted from the column in that order. Select the packing materials according to the molecular weight of the drug to be removed (e.g., Sephadex G-50 for low molecular weight, Sepharose 4B for high molecular weight).

(c) Centrifugation

When centrifugation is performed at high speeds after the addition of physiological saline (150mM NaCl), the liposomal drug is precipitated. After removal of the supernatant, add saline once again. This procedure usually needs to be repeated two or three times at 100,000rpm. Regardless of the procedure employed, it is absolutely essential to maintain both the inside and outside layers of the liposomes isotonic, and the temperature below the liposome phase transition temperature throughout the procedure. During centrifugation, the temperature should be below 40°C.

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Quantitative determination of the encapsulated drug

After removal of the non-encapsulated drug according to the methods described above, perform quantitative determination of the drug encapsulated in the liposomes after destroying the liposomes. Two methods are available for the destruction of liposomes; addition of surfactants and separation of the aqueous layer from the solvent layer by the addition of an appropriate solvent (chloroform). After decapsulating the liposomes, you can analyze the quantitative assay of the released drug by its own method.

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Nonionic Liposomes, COATSOME® EL-01-N

For enhancement of the encapsulation efficiency, the concentration of the cationic drug needs to be controlled at below 1mM. Thus, to achieve higher encapsulation efficiency, the molar ratio of DPPG (3mM) against the drug should be not less than 3.

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Improvement of the encapsulation efficiency requires maintenance of the cationic drug concentration below 5mM. If higher encapsulation efficiency is desired, the molar ratio of DPPG (15mM) against the drug should be kept at not less than 3.

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Comparison of pharmacokinetics in vivo*

| Drug Maker | Liposome | AUC(%dose/ml-min) | MRT(min) | CL tot(ml/min/Kg) |
|------------|----------------|-------------------|----------|-------------------|
| Mannitol | EL Series | 237 | 448 | 0.174 |
| Mannitol | Bangham Method | 2204 | 398 | 0.179 |
| Inulin | EL Series | 2137 | 408 | 0.175 |
| Inulin | Bangham Method | 2315 | 373 | 0.174 |

*K.Yachi et al, *Biopharm. & Drug Dispos.* **17**, 589-605(1996):Composition:DPPC:DPPG:Chol=27:53:20 (Anionic Liposomes)

Cationic Liposomes for gene delivery (COATSOME® EL-01-D)

Transfection Activities *in vitro**¹ by Luciferase Assay*⁵ (Light units/mg protein · sec)

| | HRA | MEIL | ES-2 |
|-----------|--------|-------|-------|
| Serum (-) | 6,048 | 6,291 | 1,325 |
| Serum (+) | 11,378 | 2,913 | 689 |

Transfection Activities *in vivo**² in intraperitoneal disseminated tumors*⁵

| | Percentage of LacZ-positive cells | n |
|--------------------------|-----------------------------------|---|
| COATSOME EL-01-D | 1.00 ± 0.11 | 5 |
| Commercially available A | 0.38 ± 0.26 | 3 |
| Commercially available B | 0.62 ± 0.21 | 3 |
| Commercially available C | 0.23 ± 0.23 | 3 |

Relationship Between Transfection Efficiencies and Cell Mitotic Activities*³,*⁵

| Cell lines | Percentage of LacZ-positive cells | | Labeling index (%) |
|-------------------|-----------------------------------|--------------------------|--------------------|
| | COATSOME EL-01-D | Commercially available C | |
| Cancer cell lines | | | |
| HRA | 42.9 ± 3.8 | 4.7 ± 1.0 | 68.2 ± 1.6 |
| mEiIL | 4.5 ± 1.0 | 0.6 ± 0.1 | 41.0 ± 1.2 |
| ES-2 | 23.7 ± 1.9 | 3.3 ± 0.9 | 41.9 ± 4.5 |
| OVHS-1 | 0.8 ± 0.4 | 0.9 ± 0.3 | 6.1 ± 1.4 |
| MCAS | 1.3 ± 0.1 | < 0.1 | 51.5 ± 6.1 |
| SKOV3 | 11.6 ± 1.8 | 0.2 ± 0.2 | 35.2 ± 6.7 |
| OVCAR3 | 12.3 ± 1.1 | 0.4 ± 0.4 | 34.4 ± 3.7 |
| KK 2.7 ± | 2.7 ± 0.7 | 0.1 ± 0.1 | 26.2 ± 4.3 |
| KOC-3S | 5.9 ± 2.1 | 2.6 ± 1.6 | ND* ⁴ |
| Nakajima | 4.9 ± 1.7 | 0.1 ± 0.1 | ND* ⁴ |
| KF | 15.5 ± 2.3 | 0.1 ± 0.1 | ND* ⁴ |
| SW626 | 26.8 ± 4.6 | 0.2 ± 0.1 | ND* ⁴ |
| Colo320DM | 32.3 ± 3.2 | 5.2 ± 1.1 | ND* ⁴ |
| HRA | 42.9 ± 3.8 | 4.7 ± 1.0 | ND* ⁴ |
| mEiIL | 4.5 ± 1.0 | 0.6 ± 0.1 | ND* ⁴ |
| Normal cell lines | | | |
| Fibroblast (TIG) | 9.6 ± 2.2 | 0.9 ± 0.2 | 42.4 ± 10.7 |
| Fibroblast (IMR) | 2.5 ± 0.4 | 1.6 ± 0.5 | 14.2 ± 3.2 |
| HUVEC0103 | < 0.1 | < 0.1 | 22.2 ± 1.6 |
| HUVEC1204 | < 0.1 | 0.2 ± 0.1 | 11.3 ± 0.2 |
| HUVEC0923 | < 0.1 | < 0.1 | 18.2 ± 4.1 |

*1)Transfection performed in the absence or presence of 10% fetal bovine serum (FBS). Data shown represent mean of three experiments.

*2)mEIL cells growing in peritoneal cavities of nude mice were transfected with liposome/CAG-lacZ (20μg) and the percentage of lacZ-positive cells was determined. Data shown represent mean: ± SD.

*3)Percentages of LacZ-positive cells and labeling indexes were determined. Transfection was performed in the presence of 10% FBS. Data shown represent mean ±SD. n=3

*4)ND = Not done

*5)Reference: A. Kikuchi et al., *HUMAN GENE THERAPY*, **10**: 947-955 (1999).

7. Transferrin Pre-liposome

Numerous transporters and select markers are expressed on the surfaces of cells; accordingly, with the use of liposomes chemically bound to, for example, transferrin, targeting of cancer cells over-expressing transferrin receptors might become feasible. By using the lipid mixtures available from NOF, it would be possible to prepare liposomes possessing custom cell targeting ability at your own laboratories.



TF-01-PL is a specific lipid mixture containing phospholipids conjugated with phosphatidylcholine and transferrin. If the directions in the instructions provided with your purchase are adopted, drug-containing liposomes with cancer cell targeting ability can be prepared.

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| Product name | Characteristics |
|--------------|--|
| TF-01-PL | Dried powder consists of PC and other suitable lipids containing transferrin attached phospholipid |

Transferrin-bound liposomes are directed toward the transferrin receptors on the surface of cancer cells.

- Reagent for preparing transferrin-modified liposomes

Preparation of the drug-containing liposomes involves three steps.

- Addition of the drug solution
- Mixing and hydration of the drug with the reagent
- Preparation of the liposomes by size adjustment and sterilization

Increased effectiveness of the drug

- Prolonged half-life of the drug in blood
- Avoidance of RES uptake

Facilitation of R&D of DDS preparations

- Composed of highly purified starting materials, including high-purity phospholipids
- Ready for bulk supply

This product is intended to be used for research purposes only. They are not to be used for drug or diagnostic purpose, nor are intended for human use. They shall not be used as food, cosmetics or utensils, etc.

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
Specifications:

Content: Phospholipid, cholesterol, and transferrin (Human Holo-transferrin)
Volume: Lipid: 20mg, Transferrin: 1mg per vial
Container: 10ml glass vial
Storage: Keep refrigerated

Reference:

- 1) O. Ishida et al., *Pharm. Res.* **18** (7), 1042-1048 (2001)
- 2) H. Inuma et al., *Int. J. Cancer* **99**, 130-137 (2002)

[How to prepare transferrin-modified drug-containing liposomes]

- ① Remove TF-01-PL from the refrigerator and bring to room temperature.
- ② Inject drug solution into the vial.
- ③ Gently shake the vial and mix the drug solution with the TF-01-PL powder to hydrate the powder.
- ④ Shake it gently every 5 minutes for 15 minutes (rolling shakers can also be used) so to ensure the formation of a homogeneous mixture.

Note 1) The amount of drug solution:
This product contains 20 mg lipid in a vial. Addition of the drug solution to obtain a lipid of 1-2% is recommended (the amount of drug solution to be added is 1-2 mL). Be careful not to lose the amount of solution during the filtering process described below, since the amount of solution is very small. Adding 1ml of purified water will make an approximately 30mM phosphate buffer. Adjust the drug solution buffer or salt concentration accordingly.

<Syringe filtration>

- ⑤ Suction the mixed solution prepared in ④.
- ⑥ Attach a disposable cellulose acetate filter (pore diameter: 0.22μm) to the syringe and filter the solution.
- ⑦ Transferrin-modified drug-containing liposomes (particle diameter: approximately 300nm) are obtained.

Note 2) Syringe filtration
This simple method can be employed to prepare liposomes when there are no specific requirements in terms of the liposome diameter or the amount of drug to be encapsulated in the liposomes. Aseptic injection directly into the vial after filtration will provide aseptic filtration at the same time. This product yields approximately 800μl of liposomes when 1ml of the drug solution is used and the solution is filtered with a 13-mm- diameter filter.

<High-Pressure filtration>

- ⑧ Put solution ④ into a pressure filtration system equipped with a polycarbonate filter.
- ⑨ Operate the pressure filtration system to filter the solution.
- ⑩ Transferrin-modified drug-containing liposomes are obtained. The particle diameter will be approximately 200 nm when a 0.2μm pore diameter filter is used. The drug-containing liposomes can be used after dilution or purification as needed.

Note 3) Refining the outer water phase
When the drug-containing liposomes are biologically evaluated, the drug in the outer water phase should be used after purification, according to the purpose of your research. Purification in the outer water phase can be performed by conventional methods, such as ultrafiltration, dialysis, gel filtration and centrifugation.

Note 4) Concentration
Lipid concentration of 10mg/mL (1%) or less is recommended.

Note 5) High-Pressure filtration system
Contact us for further information and purchase.

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8. Liposomal Formulations

Optimization & manufacturing of customer's liposome encapsulated hydrophilic, hydrophobic, cationic drugs

Based on our ample experience accumulated in the phospholipid/liposome fields, we are ready to comply with custom manufacture and custom development of the drug-encapsulated liposomes. We can offer proposals for the development of various functional liposomes, such as those with improved half-life of the drug in the blood and those with immunoregulatory functions.

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

| Functional Liposomes |
|--|
| Long circulation |
| Immuno-regulation (vaccine, allergy-therapy) |

| Characteristics of immuno-regulating Liposomes |
|--|
| NOF has developed immuno-regulating Liposomes |
| Inducing substantial IgG antibody |
| Inducing CTL (Cytotoxic Lymphocyte) |
| Suppressing IgE antibody |

NOF has developed immunoregulatory liposomes and technology. This technology essentially requires production of an appropriate amount of IgG antibody, while concomitantly inhibiting the production of IgE antibody. Due to the well-known advantages of strong cellular immunity (CTL activity), this technology can be applied to the preparation of vaccines for infections and cancer, as well as for the treatment of allergies, including pollen allergy.

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▲ Aluminium Adjuvant
● Liposome Adjuvant

Anti-OVA antibody production in mice with immunized Aluminum and Liposomal adjuvants. 8 weeks old female BALBC mice were immunized intraperitoneally. After immunization, blood sample were taken from the tail vein. Anti-OVA antibodies in the sera were determined.

9. LIPONIZER™ for Liposome Production


LIPONIZER™ is available only in Japan.

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


LP-142

Example of LP-142 use

10. NOF's Capabilities and Customer Advantages

GMP manufacturing facility

Bulk phospholipids produced at our own validated facility under cGMP operation, have been supplied worldwide commercial drugs for many years, with a good reputation for quality and reliability.

Drug Master Fileings

NOF has already submitted more than 32 DMFs regarding various kinds of NOF phospholipids manufactured to the US FDA, so that customers can rely on our quality. We submit DMFs according to customers' requests and development needs.

Custom Synthesis

Custom synthetic phospholipids and lipids can be produced at customer's request.

Price List for Reagents

Price list for phospholipids is downloaded from the NOF CORPORATION website: <http://www.phospholipid.jp/>. The reagents grades are high quality products as well as GMP products.

Ultra-Purity Polysorbate 80

Polysorbate 80(HX2)TM


Polyoxyethylene Sorbitan Oleate

Polysorbate 80(HX2) produced by NOF CORPORATION is a high-purity product of the highest quality in the world, with extremely low aldehyde and peroxide levels.

These characteristics contribute to its safety; in a rat study, Polysorbate 80(HX2), as compared to conventional Polysorbate 80 formulations, triggered less histamine release from rat mast cells. The latest in vivo dog study has also indicated less histamine release with our Polysorbate 80(HX2); as a natural consequence, because of the lower incidence of allergy noted, our Polysorbate 80(HX2) has drawn worldwide attention.

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
^{*1} JP:Japanese Pharmacopoeia
EP:European Pharmacopoeia
NF:National Formulary


Appearance comparison of Ultra-purity Polysorbate 80

Degranulation Test (An allergic study model)

Effect of Polysorbate 80 on degranulation of RBL-2H3 mast cells. Cells were treated with different concentrations of Polysorbate 80 for 60 mins. The degree of degranulation was determined by measurement of the released β -hexosaminidase into the supernatant.

Cell Toxicity Test


Effect of Polysorbate 80 concentration on the cytotoxicity using SIRC Cells. Cells were treated with each Polysorbate 80 for 24 hrs. The number of viable cells was determined by the Neutral Red Uptake method.

Influence on Basophilic leukocyte


Basophilic leukocytes were immersed in 1% Polysorbate 80 solution for 30 min. The cells were immobilized in glutaraldehyde and observed by SEM.

Hemolysis Test


Effect of Polysorbate 80 concentration on the Hemolytic Ratio. Red blood cells from the guinea pig were treated with Polysorbate 80 for 60 mins. The hemolytic ratio was determined by the absorbance of the solution at 576 nm.

Survival Rate

Effect of Polysorbate 80 concentration on the Survival Rate of BALB/c mice. Polysorbate 80 was diluted with PBS and administered to mice by intravenous injection.

Stability Test

This formaldehyde test was performed under N₂ sealed conditions at 40°C.

Status of European Pharmacopoeia (EP)

European Pharmacopoeia (EP) 5th edition stipulates a highest limit for oleic acid (58-85%) in Polysorbate 80. It is paradoxical that the highest-quality Polysorbate 80(HX2) in the world does not satisfy the specification of EP.

As a result of our negotiations with the EP Agency, the upper limit for the content of oleic acid has been eliminated from EP, ed. 5.4.

Now, our Polysorbate 80(HX2) meets the requirements stipulated in the Tripartite Compendia (JP, EP, NF).

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IV


Novel Excipients for Pharmaceutical Development

PUREBRIGHT® MB Series

MPC™ Polymers for Hydrophobic Drug Formulations

NOF provides unique hydrophilic polymers (PUREBRIGHT MB) that allow ready solubilization of hydrophobic drugs in aqueous solutions. These polymers consist of copolymers of 2-methacryloyloxyethyl phosphorylcholine (MPC™) and n-butyl methacrylate (BMA). Introduction of hydrophobic BMA units leads to solubilization of drugs via hydrophobic interactions.

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Advantage of PUREBRIGHT® MB Series

- High Solubility Performance for any Hydrophobic Drug
- Phosphorylcholine-Base Biocompatible Polymer
- No Irritation
- No Toxicity
- Nano Particle Drug Delivery possible

High Solubility Performance for any Hydrophobic Drug

1. High Solubility Performance

Lyophilized powder of paclitaxel with PUREBRIGHT MB-37(left) can be easily dissolved in water.

- Paclitaxel : 2mg
- PUREBRIGHT MB-37 : 30mg
- Distilled water : 1mL


Reference:

T. Konno, J. Watanabe, K. Ishihara, *J.Biomed. Mater. Res.* **65A**, 209-214 (2003)


| Product name | Description | MW |
|-----------------------|-------------|------------|
| PUREBRIGHT MB-37-50T | MPC-Polymer | ca.30,000 |
| PUREBRIGHT MB-37-100T | MPC-Polymer | ca.100,000 |

Patent No. : JP3571878, JP4233251, JP4402350, US6214957

2. Solubilization of Paclitaxel

3. Typical solubilization method

4. Relative solubilization activity of PUREBRIGHT® MB series vs. DMSO solution control

| | Abbreviation | Water-solubility |
|---------------|--------------|----------------------|
| AmphotericinB | AmB | 0.0005 mg/mL |
| Camptothecin | CTP | 0.03 mg/mL |
| Paclitaxel | TXL | Less than 0.01 mg/mL |

Reference:

1. T. Konno, et al., *J. Biomed. Mater. Res.* **65A**(2), 209-214 (2003)
2. T. Konno, et al., *Proceed. Int'l. Symp. Control. Rel. Bioact. Mater.* **29**, 464 (2002)
3. M. Wada, et al., *Anticancer Res.* **27**, 1431-1436 (2007)

| Concentration of PUREBRIGHT® MB Series | ratio of abs. from UV | | ratio of peak area from HPLC |
|---|--------------------------|------------|---------------------------------|
| | AmB | CTP | TXL |
| 10 wt% | 100.2 | 5.8 | 735 |
| 1 wt% | 84.9 | 2.8 | 618 |
| 0.1 wt% | 1.8 | 2.8 | 580 |
| 0.01 wt% | 1.4 | 1.4 | 15 |
| 10 wt% DMSO aqueous solution (Control) | 1.0 | 1.0 | 1.0 |


5. Toxicological Data

Safety data of PUREBRIGHT® MB-37 series shows similar activity in comparison with other solubilizers, such as poloxamer and polyoxyethylene castor oil. LDs of more than 25,000mg/kg (mice, orally) and more than 1,000mg/kg (rats, iv) have been obtained.

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| | |
|---------------------------|---|
| Acute toxicity (100T) | Oral LD ₅₀ (Mice) > 25,000 mg/Kg |
| Acute toxicity (50T) | Intravenous injection LD ₅₀ (Rats) > 1,000 mg/Kg |
| Antigenicity (50T) | Negative |
| Mutagenicity (50T) | Negative |
| Skin sensitization (100T) | Non skin sensitization (Guinea pig) |
| Skin toxicity (100T) | Non-primary irritation (Rabbit) |
| Eye irritation (100T) | Practically non-stimulative to eye mucosa (Rabbit) |

6. Unique Feature : Water soluble film of drug and PUREBRIGHT® MB-37 series

Reference:

1. T. Konno, et al., *J. Biomed. Mater. Res.* **65A**(2), 209-214 (2003)
2. T. Konno, et al., *Proceed. Int'l. Symp. Control. Rel. Bioact. Mater.* **29**, 464 (2002)

7. Pharmaceutical applications (published)

1. Oncology Drug Formulation
 - 1) H. Takeuchi et al., *Digestion* **82**, 187-191 (2010)
 - 2) T. Kamei et al., *Cancer Sci.* **102**(1), 200-205 (2010)
 - 3) D. Soma et al., *Cancer Sci.* **100**(10), 1979-1985 (2009)
2. DNA Delivery
 - M. Ukawa et al., *Biomaterials* **31**(24), 6355-6362 (2010)


Self-emulsifying Solubilizer SL-11

Novel Nano Formulation for Hydrophobic Drugs

Solubilizer SL-11, occurring as a transparent liquid, can be used to prepare nano-emulsions of various hydrophobic drugs (particle sizes are below 50 nm). Since this possesses excellent self-emulsifying properties, this is the most appropriate agent for SEDDS (Self Emulsifying Drug Delivery Systems). Following oral administration, this solubilizer SL-11 augments intestinal absorption of drugs. Once the hydrophobic drug is dissolved in SL-11 as directed in the attached protocol, SEDDS reagents containing the drug can be easily prepared. NOF has also capability for custom development by using this technology.

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| Product Name | Content |
|-------------------|--------------------------------------|
| Solubilizer SL-11 | Mixture of pharmaceutical excipients |

Remarks: *Temperature may be changed depending upon drug stability.

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Cutting-Edge Technology for Self-emulsifying Agents (Particle size: <50nm)

- Easy preparation
- Applicable to various kinds of hydrophobic drugs
- Applicable to soft capsules for oral use
- Most suitable for SEDDS agent
- Improvement of intestinal absorption following oral administration
- Custom development for formulations
- Comply with bulk supply requests

[Specification]

Content : mixture of pharmaceutical excipients
Volume : 50g , 100g , 1kg
Container : glass bottle
Storage : Storage under room temperature.
Keep container closed.

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
PUREBRIGHT® SL Series

Hydrophobic Drug Solubilization Kit

NOF produces Solubilization Kits for hydrophobic drugs using our various DDS technologies. The procedures for the use of these kits are very simple, facilitating your screening of new chemical entities at early test stages.

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PUREBRIGHT®


| Product name | Description |
|-------------------------------------|--|
| PUREBRIGHT SL-110 ¹⁾⁻²⁾ | MPC™ polymer based (PUREBRIGHT MB-37-100T) |
| PUREBRIGHT SL-220 | PEG-phospholipid based (SUNBRIGHT DSPE-020CN) |
| PUREBRIGHT SL-310 ³⁾⁻¹⁰⁾ | PEG-cholesterol based1 (SUNBRIGHT CS-010) |
| PUREBRIGHT SL-350 ³⁾⁻¹⁰⁾ | PEG-cholesterol based2 (SUNBRIGHT CS-050) |
| PUREBRIGHT SL-411 | Phospholipid complex based (COATSOME MC-4040/COATSOME MG-4040LS) |

SUNBRIGHT® COATSOME® PUREBRIGHT®

References:

- 1) K. Ishihara et al., *Polym.J.* **31**, 1231-1236 (1999)
- 2) T. Konno et al., *J.Biomed.Mater.Res.* **65A(2)**, 209-214 (2003)
- 3) H. Ishiwata et al., *Chem. Pharm. Bull.* **46**, 1907-1913 (1998)
- 4) Uchegbu IF. et al., *Int. J. Pharm.* **155**, 7-17 (1997)
- 5) Tasset C. et al., *Int. J. Pharm.* **58**, 41-48 (1990)
- 6) Uchegbu IF. et al., *Pharm. Res.* **12**, 1019-1024 (1995)
- 7) Uchegbu IF. et al., *S.T.P. Pharma Sci.* **6**, 33-43 (1996)
- 8) Dimitrijevic D. et al., *J. Pharm. Pharmacol.* **49**, 611-616 (1997)
- 9) Uchegbu IF. et al., *J. Pharm. Pharmacol.* **49**, 606-610 (1997)
- 10) Auth R. et al., *Akt. Dermatol.* **10(6)**, 215-220 (1984)

Hydrophobic Solubilization Methods

< Preparation >

The following instructions are for preparation of 5 vials of the test solution simultaneously.

(1) Weigh drugs into test tube or cuvette according to the desired concentration in the final solution. If you wish to make 1mg/mL of the final solution, weigh 12 mg of the drug, add 3mL of ethanol (select other solvent when drug is not soluble in ethanol) into the test tube, and then dissolve completely by gentle shaking.

(2) Pipette out 500μL of the drug solution, and transfer the solution into each vial of PUREBRIGHT SL Kit respectively.

(3) Dissolve completely by shaking. Warm the vial to 50°C* to obtain a clear solution.

(4) Evaporate the solvent from the vial according to one of the following methods

a) Place the vial into vacuum chamber previously warmed to 50°C*, then evaporate the solvent under vacuum.

b) Evaporate the solvent by a N₂ gas stream at 50°C* over vials placed in a warm bath. Exercise caution in handling to avoid fires or accidents caused by volatile solvent.

(5) Add 2mL of deionized water into the vial, and the dissolve completely with shaking.

Add 1.8mL of deionized water in case of using PUREBRIGHT SL-411.

Drug in aqueous solution (e.g., 1mg/ml) is now ready.

Remarks: *Temperature may be changed depending upon the drug stability.

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Drug Solubilization with PUREBRIGHT® (薬剤の可溶化例)

| Drug | Pep | IND | MPS | CLF | NFD | PRB | CLM | PTX |
|------------------------|-----|-----|------|-----|-----|-----|-----|-----|
| Concentration (mg/mL) | 6.0 | 2.5 | 0.75 | 7.5 | 2.5 | 5.0 | 0.5 | 1.0 |
| PUREBRIGHT SL-110 | I | I | S | I | I | A | I | S |
| PUREBRIGHT SL-220 | S | S | N | I | I | I | I | S |
| PUREBRIGHT SL-310 | S | S | S | S | S | I | S | S |
| PUREBRIGHT SL-350 | I | I | N | N | I | I | I | S |
| PUREBRIGHT SL-411 | S | N | N | N | N | S | N | N |
| Polyoxyl 35 Castor oil | I | I | I | A | I | I | I | I |

PUREBRIGHT & Polyoxyl 35 Castor oil : 3% (w/w) S: Soluble A: Almost soluble I: Insoluble N: No data

Pep: Peptide type Drug

IND: Indomethacin

MPS: Medroxyprogesterone acetate

CLF: Clofibrate

NFD: Nifedipine

PRB: Probucol

CLM: Chlormadinone acetate

PTX: Paclitaxel


NOFABLE™ Series

Ultra-Purity Oleic Acid and their Derivatives


NOF has established the manufacturing processes for 99% oleic acid purity in commercial production scale. Removal of unsaturated fatty acids susceptible to oxidation is eventually associated with improvement of product stability against oxidation, with the result that derivatives with less impurities (peroxides and aldehydes) can be obtained. Ultra-pure oleic acid derivatives are the most suitable for efficacious drug formulations.

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

Gas chromatogram for NOF Oleic acid compared to conventional product


Crystal of EXTRA OLEIN 99

Oxidation Stability of Oleic Acid

Formation of peroxides in oleic acid during incubation at 50°C.


Dermal Absorption Enhancement with Oleic Acid

Dermal Absorption Enhancement of Indomethacin with Highly Purified Oleic Acid (EXTRA Series).

2. Ultra-Purity Oleic Acid Derivatives

•Sorbitan Oleate


| Product name | Description | Oleic acid Content | CAS NO. | Regulatory Status*1 |
|-----------------|---------------------|--------------------|-----------|---------------------|
| NOFABLE™ SO-991 | Sorbitan monooleate | >99% | 1338-43-8 | JPE |

•Glycerol Oleate

| Product name | Description | Oleic acid Content | CAS NO. | Regulatory Status*1 |
|-----------------|---------------------|--------------------|------------|---------------------|
| NOFABLE™ GO-991 | Glycerol monooleate | >99% | 25496-72-4 | JPE |

•Ethyl Oleate

| Product name | Description | Oleic acid Content | CAS NO. | Regulatory Status*1 |
|----------------|--------------|--------------------|----------|---------------------|
| NOFABLE™ EO-99 | Ethyl oleate | >99% | 111-62-6 | JPE |

*1 JPE : Japanese Pharmaceutical Excipients

SUNBRIGHT® DKH-02HB, DKH-03HB and DKH-04HB

MACROGOL (PEG-200, 300 and 400)

NOF produces top-quality Polyethylene Glycols containing only with extremely low levels of EG (ethylene glycols) and DEG (diethylene glycols), whereby our Polyethylene Glycols can be specifically used as pharmaceutical excipients. From the toxicological perspective, EP has stipulated that the upper limits of EG and DEG should not exceed 4000 ppm, while the USP stipulates that the upper limits of EG and DEG should not exceed 2500 ppm. NOF can guarantee much lower levels of EG and DEG than those stipulated by both EP and USP. In addition, we pride ourselves in the fact that our Polyethylene Glycols contain extremely low levels of impurities, including ethylene oxides, dioxane, peroxides and aldehyde derivatives. NOF can supply high-quality PEG 200, 300 and 400, all of which would meet the specifications of EP, USP, JP and JPE.

Concerning this product, NOF has officially registered its product patent in the United States of America.

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

| | Product Name | Description | MW |
|--------|--------------------|--------------------|-----|
| PEG200 | SUNBRIGHT DKH-02HB | Polyethyleneglycol | 200 |
| PEG300 | SUNBRIGHT DKH-03HB | Polyethyleneglycol | 300 |
| PEG400 | SUNBRIGHT DKH-04HB | Polyethyleneglycol | 400 |

SUNBRIGHT®

Stability Test


This formaldehyde test was performed under N₂ sealed condition at 40°C.

This pH test was performed under N₂ sealed condition at 40°C.

EG/DEG Content

Specifications

| ITEM | Method | Specification | | |
|--|--------|--|----------|----------|
| | | DKH-02HB | DKH-03HB | DKH-04HB |
| Completeness and color solution | | Pass | | |
| Viscosity | | 3.9-4.8 | 5.4-6.4 | 6.8-8.0 |
| Average molecular weight | | 190-210 | 285-315 | 380-420 |
| pH | | 4.5-7.5 | | |
| Residue on ignition (%) | USP | NMT 0.1 | | |
| Heavy metals (ppm) | | NMT 5 | | |
| Free ethylene oxide and 1,4-dioxane (ppm) | | NMT 10 (Ethylene oxide : NMT 1, 1,4-Dioxane : NMT 10) | | |
| Limit of ethylene glycol and diethylene glycol (%) | | NMT 0.05 | | |
| Organic volatile impurities | | Pass | | |
| Peroxide value (meq/Kg) | NOF | NMT 3 | | |
| Formaldehyde (ppm) | | NMT 5 | | |


NOF can also supply products without anti-oxidant agent, BHT.

SUNBRIGHT® OE Series

Biocompatible PEG Anchors for Cell Membrane Insertion [BAM]

NOF supplies a unique Biocompatible Anchor for cell Membranes (BAM), without eliciting observable cell damage. BAM contains both oleyl groups as a hydrophobic cell membrane anchor and polyethylene glycol (PEG), which imparts hydrophilicity. Its chemical structure is designed to have various reactive groups at the PEG terminals allowing them to attach to physiologically active substances or surfaces of materials. Employment of BAM enables researchers to modify the surfaces of cells and tissues with physiologically active substances, such as proteins or drugs without injuring cells and tissues, and also to immobilize live cells and tissues on the surfaces of various kinds of materials. Since BAM exerts excellent cell interfacing functions beyond conventional concepts, it can be extensively used in pharmaceuticals and cosmetics.

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
| Product name | Reactive Group (X) | Mw |
|--------------------|--|-------|
| SUNBRIGHT OE-020CS | | 2,000 |
| SUNBRIGHT OE-040CS | -CO-CH ₂ CH ₂ -COO-NHS | 4,000 |
| SUNBRIGHT OE-080CS | | 8,000 |

SUNBRIGHT®


Application I :

Anchoring of BAM-Protein conjugate to the cell membrane

Application II :

Anchoring of cell on a BAM modified surface

References:

- 1) K. Kato, et. al, *Biotechnol. Prog.* **20**, 897-904 (2004)
- 2) K. Kato, et. al, *BioTechniques* **35**,1014-1021 (2003)

Ordering Information

| | |
|---|--|
| In the Americas <u>NOF AMERICA CORPORATION</u> One North Broadway, Suite 1012 White Plains, NY 10601, USA Tel : 914-681-9790 Fax : 914-681-9791 E-mail : info@nofamerica.com Home Page : http://www.nofamerica.com | In the region except in Americas, Europe <u>NOF CORPORATION</u> DDS Development Division 20-3 Ebisu 4-chome, Shibuya-ku Tokyo, 150-6019, Japan Tel : +81 3 5424 6741 Fax : +81 3 5468 0901 E -mail : ddsinfo@nof.co.jp Home Page : http://www.dds-drug.com |
| In Europe <u>NOF EUROPE (BELGIUM) N.V.</u> Bouwelven 1 BE-2280 Grobbendonk, Belgium Tel : +32 14 25 98 26 Fax : +32 14 22 45 63 E -mail : info@nofeurope.com Home Page : http://www.nofeurope.com | In Japan (日本国内) <u>日油株式会社</u> DDS 事業部 〒150-6019 東京都渋谷区恵比寿 4-20-3 恵比寿ガーデンプレイスタワー Tel : 03-5424-6741 Fax : 03-5468-0901 E -mail : ddsinfo@nof.co.jp Home Page : http://www.dds-drug.com |
| A faxed confirmation is required for every order | |

Delivery : Standard deliveries in Japan in-stock product are by domestic courier. International shipment is delivered via air transport. The buyer is responsible for all international duties, tariffs, taxes and other importing fee. Except in a few case, we have a stock of the products in this catalogue without made-to-order product indicated with asterisk. In case we don't have any stock of the product, the lead time of its shipment will be 6 to 8 weeks after buyer's purchase order.

Warranty: Unless otherwise any consent between NOF and buyer, all of the products in this catalogue are offered for research purpose only and they are not intended for human use. The buyer assumes all risk and liability for the use and/or results obtained by the use of products described in this catalogue.

Custom Orders : Order of derivatives of other molecular weight and/or chemical structure like PEG-PLA and so on require custom synthesis. Custom orders are not returnable or refundable. Please contact our sale office above shown for a quote on your request product.

Pricing Offer : Pricing is available for orders of 1 gram or more for non-GMP. GMP pricing is available for orders of 1000 grams or more. Please call our sale office above shown for each quote.

GMP manufacturing : Many of the products listed in this catalogue are currently available as non-GMP materials. However, some products are already manufactured to cGMP at dedicated cGMP plant. Please ask our sale office for the detail information.

Terms and Conditions of Sales

Definitions: “Seller” means NOF CORPORATION; “Buyer” means the purchaser of Products hereunder; and “Products” means the items or goods sold by Seller to Buyer hereunder.

Acceptance: All sales are subject and limited to, and conditioned upon, these Term and Conditions. No variations bind Seller unless contained in a writing signed by an offer of Seller.

Changes: Orders may be changed or amended only pursuant to a writing signed by both Buyer and Seller setting forth the changes to be made. Buyer may not cancel an order unless Seller agrees to such cancellation in writing. In such event, Buyer shall pay all storage and shipment costs, costs of producing non-standard material, costs of purchasing non-returnable materials, cancellation costs imposed on Seller by its suppliers, and any other costs resulting from such cancellation, all as specified to Buyer by Seller.

Delivery: All sales are Delivered Duty Unpaid (DDU, International Rules for the Interpretation of the Trade Terms 2000) at the port of destination Seller shall deliver the Products to Buyer at Seller's cost and Seller bears all risk of loss or damage in transit thereafter Seller reserves the right to make delivery in installments, all such installments shall be separately invoiced and paid for when due per invoice, without regard to subsequent deliveries. Delay in Delivery of any installment shall not relieve Buyer of its obligation to accept remaining deliveries.

Inspection: Within five (5) days after receipt, (a) Buyer shall inspect whether Products contains for chipped glass, broken seals and other patent defects, and (b) shall notify Seller in writing of any claims for shortages, defects or damage identified in such inspection. Thereafter, Buyer shall hold Products for Seller pursuant to written instructions concerning disposition. If Buyer fails to notify Seller within such five (5) day period, Products shall conclusively be deemed to conform to the terms and conditions of the sales and transactions covered herein (hereinafter “Sales”) with respect to any such shortages, defects or damage and to have been irrevocably accepted by the Buyer in their short, defective or damaged condition.

Force Majeure: Seller shall not be liable for any delay or other failure to perform hereunder due to any cause beyond Seller's reasonable control, including, without limitation, acts of Buyer, governmental actions or regulation, fire, explosion, accident, theft, vandalism, riot, acts of war, strikes or other labor difficulties, lighting, flood, tornado, typhoon, windstorm or other acts of God, transportation delays, or inability to obtain necessary labor, fuel, materials, supplies or power at current prices.

Allocation of Products: If Seller is unable for any reason to supply the total demands for Products specified in Buyer's order, Seller may allocate its available supply among any or all buyers on such basis as Seller may deem fair and practical, without liability for any failure of performance which may result therefrom.

Payment: Terms of payment are net 30 days of date of invoice. If insecure as to the collection of amounts due, Seller may delay or postpone delivery of Products and change the payment terms to payment in full or in part in advance of shipment of the undelivered Products. If Buyer defaults in the payment of the purchase price or other terms of this or any other order, Seller may defer delivery, cancel the Sales, or sell products on hand for the Buyer's account and apply the proceeds against the purchase price. Buyer shall pay any balance to Seller on demand. Buyer shall pay all costs, including, without limitation, reasonable attorney and accounting fees and other expenses of collection resulting from any default by Buyer in any of the terms hereof. A late payment charge of 1.5% per month shall accrue on delinquent amounts until paid.

Taxes: Buyer shall pay any tax, duty, custom, fee or charge imposed by any governmental authority on or measured by the Sales, and Buyer shall reimburse Seller for any payment by Seller thereof.

Pricing: Buyer shall remit payment by wire transfer to the bank account mentioned in invoice. Buyer shall contact Seller for current prices if Buyer requires this information prior to placing an order. Seller guarantees pricing contained in its written quotations for 30 days. When placing an order, Buyer should reference Seller's quoted prices or the quote number. Seller will contact Buyer by phone, fax or e-mail in the event of any discrepancy between Buyer's order and Seller's pricing or other terms.

Warranties: Seller will make available to Buyer a Certificate of Analysis for Products purchased. Provided Buyer notifies Seller in writing of any shortages or patent defects or damage within the five (5) day period described in the “Inspection” section above, or, for any other products defect, within thirty (30) days after the date of the invoice for such Products, Seller will reprocess or replace, at Seller's sole option, any Products which are defective in material or workmanship. THE FOREGOING CONSTITUTES BUYER'S SOLE RECOURSE. SELLER EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES AND LIABILITY WHATSOEVER, WHETHER EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY AND/OR FITNESS FOR A PARTICULAR PURPOSE, BUYER BEARS ALL RISK RESULTING FROM THE PARTICULAR USE IT CHOOSES FOR THE PRODUCTS. No employee, agent or representative of Buyer has the authority to bind Seller to any oral representation or warranty concerning any Products. Seller's warranties shall not be effective if Seller determines that Buyer has misused Products in any manner, failed to use Products in accordance with industry standards and practices, or failed to use Products in accordance with instructions, if any, furnished by Seller.

Seller's sole and exclusive liability and Buyer's exclusive remedy with respect to Products proved to Seller's satisfaction to be defective or nonconforming shall be the replacement or reprocessing of such Products without charge or refund of the purchase price, in Seller's sole discretion, upon the return of such Products in accordance with Seller's instructions. SELLER SHALL NOT IN ANY EVENT BE LIABLE FOR INDIRECT, INCIDENTAL, CONSEQUENTIAL, PUNITIVE OR SPECIAL DAMAGES OF ANY KIND RESULTING FROM THE SALE, DELIVERY, USE OR FAILURE OF THE PRODUCTS, EVEN IF SELLER HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGE INCLUDING, WITHOUT LIMITATION, LIABILITY FOR LOSS OF USE, LOSS OF WORK IN PROGRESS, DOWN TIME, LOSS OF REVENUE OR PROFITS, FAILURE TO REALIZE SAVINGS, LOSS OF PRODUCTS OF BUYER OR OTHER USE OR ANY LIABILITY OF BUYER TO A THIRD PARTY ON ACCOUNT OF SUCH LOSS, OR FOR ANY LABOR OR ANY OTHER EXPENSE, DAMAGE OR LOSS OCCASIONED BY SUCH PRODUCTS INCLUDING PERSONAL INJURY OR PROPERTY DAMAGE UNLESS SUCH PERSONAL INJURY OR PROPERTY DAMAGE IS CAUSED BY SELLER'S GROSS NEGLIGENCE. SELLER'S AGGREGATE LIABILITY FOR DAMAGES HEREUNDER, WHETHER IN CONTRACT, TORT, OR OTHERWISE, SHALL NOT EXCEED THE PURCHASE PRICE PAID BY BUYER FOR THE PARTICULAR PRODUCTS INVOLVED. All claims must be brought within one (1) year of delivery, regardless of their nature.

Compliance with Laws: Seller certifies that to the best of its knowledge Products are produced in compliance with applicable requirements of the Labor Standard Law in Japan, as amended, and regulations, rules and orders issued pursuant thereto.

Buyers' Representations and Indemnity: Buyer represents and warrants that it shall use all Products purchased herein in accordance with “Buyer's Use of Product” as described below, and that any such use of Products will not violate any law or regulation. Buyer shall indemnify and hold harmless Seller, its employees, agents, successors, officers, and assigns, from and against any suites, losses, claims, demands, liabilities, costs and expenses (including attorney and accounting fees) that Seller may sustain or incur as a result of any claim against Seller based upon negligence, breach of warranty, strict liability in tort, contract, or any other theory of law brought by Buyer, its officers, agents, employees, successors or assigns, by Buyer's customers, by end users, by auxiliary personnel (such as freight handlers, etc.) or by other third parties, arising out of, directly or indirectly, the use of Products, or by reason of Buyer's failure to perform its obligations contained herein. Buyer shall notify Seller in writing within fifteen (15) days of Buyer's receipt of knowledge of any accident or incident involving Products which results in claimed personal injury or damage to property, and Buyer shall fully cooperate with Seller in the investigation and determination of the cause of such accident and shall make available to Seller all statements, reports and tests made by Buyer or made available to Buyer by others. The furnishing of such information to Seller and any investigation by Seller of such information or incident report shall not in any way constitute any assumption of any liability for such accident or incident by Seller.

Confidentiality: The Buyer shall keep confidential and shall not publish or otherwise disclose and shall not use for any purpose other than the purpose of the purchase of Products all confidential or proprietary information, data, documents or other materials supplied by the Seller in connection with the Sales and marked or otherwise identified in written form as “Confidential” including information derived from a site visit to Seller's facility by Buyer (hereinafter “Confidential Information”). The Buyer shall use at least the same standard of care as it uses to protect its own confidential information to ensure that the Buyer (i) only makes use of Confidential Information for purposes of the purchase of Products and (ii) does not disclose or make any unauthorized use of such Confidential Information. Notwithstanding the foregoing, the provisions of these Terms and Conditions shall not apply to Confidential Information that the Buyer can show by written records: (a) is in the public domain other than by acts of the Buyer in contravention of these Terms and Conditions; (b) was permitted to be disclosed by prior consent of the Seller; (c) has become known to the Buyer by a third party, provided such Confidential Information was not obtained by such third party directly or indirectly from the Seller under these Terms and Conditions on a confidential basis; (d) prior to disclosure, was already in the possession of the Buyer; (e) is independently developed without use of Confidential Information; and (f) is required to be disclosed to comply with any applicable law, regulation or court order, provided that the Buyer shall provide prior written notice of such disclosure to the Seller and take reasonable and lawful actions to avoid or minimize the degree of such disclosure.

Technical Assistance: At Buyer's request, Seller may, at Seller's discretion, furnish technical assistance and information with respect to Products. SELLER MAKES NO WARRANTIES OF ANY KIND OR NATURE, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE, WITH RESPECT TO TECHNICAL ASSISTANCE OR INFORMATION PROVIDED BY SELLER OR SELLER'S PERSONNEL. ANY SUGGESTIONS BY SELLER REGARDING USE, SELECTION, APPLICATION OR SUITABILITY OF THE PRODUCTS SHALL NOT BE CONSTRUED AS AN EXPRESS OR IMPLIED WARRANTY UNLESS SPECIFICALLY DESIGNATED AS SUCH IN A WRITING SIGNED BY AN OFFICER OR OTHER AUTHORIZED REPRESENTATIVE OF SELLER.

Buyer's Use of Products: Buyer agrees and acknowledges that Seller's products are intended primarily for laboratory research purpose and, unless otherwise states on product labels, Sellers's catalogue or in other literature furnished to Buyer, are not to be used for any other purposes, including, without limitation, in vitro diagnostic purposes, in foods, drugs, medical devices or cosmetics for humans or animals or for commercial purposes. Seller has not tested the products for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use, unless otherwise stated in Seller's literature furnished to Buyer.

Infringement: The Buyer shall indemnify and hold the Seller harmless from and against any and all infringement with regard to including but not limited to patent, utility model, design, trademark, copyright, trade secret or any other intellectual property right in connection with the Products in any country (hereinafter referred to as the “Infringement”). The Buyer shall, at its own cost, defend, settle or otherwise dispose of claims, suits or legal proceeding in connection with Infringement.

Returns: Products may not be returned for credit except with Seller's permission, and then only in strict compliance with Seller's return shipment instructions.

Miscellaneous: Seller's waiver of any default by Buyer, or any failure by Seller to strictly enforce any provisions of these Term and Condition or to exercise any right arising hereunder, shall not constitute a waiver of (a) any subsequent default, or (b) Seller's right to strictly enforce such Terms and Conditions or to exercise such right thereafter. All rights and remedies under the Sales are cumulative and are in addition to any other rights and remedies Seller may have at law or in equity. If any provision of these Terms and Conditions shall be held to be invalid, illegal or unenforceable, the validity, legality and enforceability of the remaining provisions shall not be affected or impaired thereby. The paragraph headings herein are for convenience only; they form no part of these Terms and Conditions and shall not affect their interpretation. These Terms and Conditions shall be binding upon, inure to the benefit of, and enforceable by, the parties hereto, and their respective heirs, personal representatives, successors and assigns.

Governing Law and Jurisdiction: All disputes as to the legality, interpretation, application, or performance of the Sales or any of these Terms and Conditions shall be governed by the law of Japan and subject to the exclusive jurisdiction of the Tokyo District Court for first instance.

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- SUNBRIGHT, COATSOME and PUREBRIGHT are registered trademarks of NOF CORPORATION in the U.S. and registered trademarks or trademarks in other countries.
- NOFABLE is a trademark of NOF CORPORATION in the U.S. and a registered trademark or trademark in other countries.
- SAINT is a trademark of Synvolux Therapeutics B.V. in the U.S. and other countries.
- LIPONIZER is a trademark of Nomura Micro Science Co.,Ltd. in the U.S. and a registered trademark or trademark in other countries.
- CIMZIA is a registered trademark of USB Group of Companies in the U.S. and a registered trademark or trademark in other countries.
- PEGASYS is a registered trademark of Hoffman-La Roche Ltd. in the U.S. and a registered trademark or trademark in other countries.
- MACUGEN is a registered trademark of Eyetech Inc. in the U.S. and a registered trademark or trademark in other countries.
- ALL other registered trademarks and trademarks are the property of their respective owners.

NOF CORPORATION

DDS Development Div.

Head Office

Yebisu Garden Place Tower, 20-3 Ebisu 4-chome, Shibuya-ku, Tokyo 150-6019

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